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DEVICE DESCRIPTION

Amplatz Renal Dilation Set-S contain:

- 8Fr guiding catheter
- 10Fr radiopaque dilation sheath
- Radiopaque Fascial dilators 6, 8 & 10 Fr of length 20 cm
- Renal dilators from 12, 14, 16, 18, 20, 22, 24, 26, 28 & 30 Fr of length 30 cm
- Teflon Amplatz sheath for dilators 24, 26, 28 & 30 Fr of length 17 cm.

Amplatz Renal Dilation Set – F contain:

- Scalpel Blade
- Trocar Needle with sheath 6 Fr
- 8 Fr guiding Catheter
- 0.035 inch PTFE extra stiff J tip guidewire
- 10 Fr Radiopaque Dilation Sheath
- Radiopaque Fascial dilators 6, 8 & 10 Fr of length 20 cm
- Renal dilators from 12, 14, 16, 18, 20, 22, 24, 26, 28 & 30 Fr of length 30 cm
- Teflon Amplatz sheath for dilators 24, 26, 28 & 30 Fr of length 17 cm.

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- Scalpel Blade
- Trocar Needle with sheath 6 Fr
- 8 Fr Guiding Catheter
- 0.035” PTFE extra stiff J tip guidewire
- 10 Fr Radiopaque Dilation Sheath
- Radiopaque Fascial dilators 6, 8 & 10 Fr of length 20 cm
- Radiopaque Renal dilators from 12, 14, 16, 18, 20 & 22 Fr of length 30 cm
- Teflon Amplatz sheath for dilators 14, 16, 18, 20 & 22 Fr of length 17 cm

INTENDED PURPOSE & DURATION

Amplatz Renal Dilation Set is intended for stretching or enlarging the fascial tissue prior to a nephrostomy procedure and to introduce the Amplatz sheath which is intended for percutaneous dilation of tissue surrounding a dilator previously placed. It is used as a working channel to maintain the established nephrostomy tract for introduction of instruments and catheters during urological procedures up to a period of 60 minutes.

INDICATIONS FOR USE & CONTRAINDICATIONS

This device is contraindicated in the presence of conditions which create unacceptable risk during percutaneous nephrostomy.

POTENTIAL COMPLICATIONS

Not limited to Tissue trauma, Bleeding, Perforation, Infection, pain.

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WARNINGS AND PRECAUTIONS

- This device is intended to be used under imaging guidance only.
- Use with caution when advancing the sheath and dilators to avoid tissues and orifice damages. Use only compatible advancing sheath and wireguide.
- Do not use if product is damaged or expired, reuse and resterilization of this device leads to loss in structural integrity and performance.
- This device is intended for use by trained physicians and this device should not be used without comprehensive knowledge of the indications, techniques and risks of the procedure.
- On completion of the procedure, dispose off the device as per local regulatory guidelines for bio-hazardous medical waste disposal.

PROCEDURE

1. Perform initial access with puncture needle, place a guidewire and perform sequential dilation as appropriate using standard procedure.
2. Perform initial access with puncture needles and place a guidewire for the purpose of sequential dilation using standard procedure.
3. Perform initial dilation using fascial dilators sequentially from smallest to the largest dilator.
4. After that insert the 10Fr dilatation sheath over the guidewire.
5. Load the 8Fr guiding catheter on the previously loaded 10Fr dilatation sheath and remove the sheath.
6. After using the guiding catheter, perform sequential dilation using renal dilators from the smallest to the largest renal dilator as per the requirement of the procedure.
7. When dilating the tract, the 8 Fr stylet must be held in position and the dilators passed into the tract using a twisting motion.
8. Now introduce Amplatz sheath of required size over the last dilator used.
9. Remove the dilator, leaving the amplatz sheath in position.
10. This can now be used as a working channel to perform the procedure.
11. Remove the device carefully on completion of the procedure.
12. Use package supplied accessories only.

SUPPLY AND STORAGE

The device is packed with primary pack made of polyethylene and the secondary package is made of polyethylene and Tyvek pouch. Supplied sterilized by Ethylene Oxide gas in peel open package.

Device shall be stored at temperature range from 20^oC to 40^oC and away from direct sunlight and other adverse environmental conditions.

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
Non-pyrogenic

Apirógeno
Apirogeno

Single sterile barrier system with protective package inside

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Blueneem, the Blueneem Logo are Registered and Trademarks of Blue Neem Medical Devices Private Limited Company.

© 2024 All rights reserved

BLUE NEEM MEDICAL DEVICES PRIVATE LIMITED

Plot Nos 270 & 271, Road No 5,
Harohalli Industrial Area, II Phase,
Kanakapura Taluk, Ramanagara -
562112, India.

Phone: +91 9739972855 / 54

E-Mail: contact@blueneem.com

Website: www.blueneem.com

**Authorized Representative in the
European Community:**

Arazy Group (Ireland) Ltd

19 Baggot Street Lower,
Dublin 2, D02 X658, Ireland

E-mail: Ireland@arazygroup.com