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Instructions for Use

EN

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Instructions for use

ES

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IT

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DEVICE DESCRIPTION

Amplatz Sheath a slender tube like device made of medical grade plastic material with tapered bevel tip at the distal end.

INTENDED PURPOSE & DURATION

Intended for percutaneous dilation of tissue surrounding a dilator previously placed. It is used as a working channel to maintain a previously established nephrostomy tract for introduction of instruments and catheters during urological procedures up to a period of 60 minutes.

INDICATIONS FOR USE & CONTRAINDICATIONS

None known

POTENTIAL COMPLICATIONS

Not limited to bleeding and renal puncture.

WARNINGS AND PRECAUTIONS

- Use with caution when advancing the sheath and dilators to avoid tissues and orifice damages. Use only compatible advancing sheath.
- Do not use if product is damaged or expired, reuse and resterilization of this device leads to loss in structural integrity and performance.
- This device is intended for use by trained physicians and this device should not be used without comprehensive knowledge of the indications, techniques and risks of the procedure.
- On completion of the procedure, dispose off the device as per local regulatory guidelines for bio-hazardous medical waste disposal.

PROCEDURE

1. After dilating the tract with chosen size of dilator, advance the amplatz sheath over the dilator slowly with firm twisting motion, until it enters into the renal pelvis using imaging guidance.
2. Once the position is confirmed remove the dilators leaving the sheath in place, for establishing working tract.
3. On completion of the procedure, gently withdraw the device.

SUPPLY AND STORAGE

The device is packed with primary pack made of polyethylene and the secondary package is made of polyethylene and Tyvek pouch. Supplied sterilized by Ethylene Oxide gas in peel open package.

Device shall be stored at temperature range from 20⁰C to 40⁰C and away from direct sunlight and other adverse environmental conditions.

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
Non-pyrogenic

Apirógeno
Apirogeno

Single sterile barrier system with protective package inside

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Blueneem, the Blueneem Logo are Registered and Trademarks of Blue Neem Medical Devices Private Limited Company.

© 2024 All rights reserved

BLUE NEEM MEDICAL DEVICES PRIVATE LIMITED

Plot Nos 270 & 271, Road No 5,
Harohalli Industrial Area, II Phase,
Kanakapura Taluk, Ramanagara -
562112, India.

Phone: +91 9739972855 / 54

E-Mail: contact@blueneem.com

Website: www.blueneem.com

**Authorized Representative in the
European Community:**

Arazy Group (Ireland) Ltd

19 Baggot Street Lower,
Dublin 2, D02 X658, Ireland

E-mail: Ireland@arazygroup.com