Bluthane[®] Ureteral Stents and Stent Sets

Instructions for Use

EN

Bluthane[®] Ureteral Stents and Stent Sets

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Bluthane® Ureteral Stents and Stent Sets

DEVICE DESCRIPTION

Bluthane® Ureteral stents are flexible tubular devices made from radiopaque polyurethane. One or both ends of the ureteral stent consist of curled pigtail self retaining loops. Both ends of the stent are open or one end of stent is closed based on the type chosen. Ports are provided along the length of the stent including the pigtail loops of the stent. The stent is marked along the length to enable visualization, advancement and placement.


If the Stent provided with monofilament string which helps in device positioning during insertion and also if stent provided with surface hydrophilic coat, this coating when activated will make stent outer surface slippery for easy insertion.

Bluthane® Ureteral Stent package includes:

- Radiopaque double pigtail stent
- Pusher/Stent Positioner with radiopaque tip
- Clamp

Bluthane® Ureteral Stent Set package includes(set components may vary)

- Radiopaque double pigtail stent
- Pusher/Stent Positioner with or without radiopaque tip
- Guidewire
- Clamp

INTENDED PURPOSE & DURATION

Bluthane® Ureteral stents are non-bioabsorbable medical device intended for treatment of obstructed ureter for drainage of urine from renal pelvis to urinary bladder and is intended for long-term use up to 180 Days.

INDICATIONS FOR USE & CONTRAINDICATIONS

Indicated to be administered in an obstructed ureter (e.g., from stone, stricture, compression, fibrosis, trauma).

Not to be implanted in conditions of extreme strictures caused by tumors.

POTENTIAL COMPLICATIONS

Potential complications include but not limited to Ureteral reflux, Urinary tract infection (Uro sepsis), Hematuria, Edema, Migration, Dysuria, Encrustation and stent fragmentation.

Bluthane® Ureteral Stents and Stent Sets

WARNINGS AND PRECAUTIONS

- Do not reuse the device.
- Do not use if package is opened or damaged. Do not use the product if there is suspect that the device is not sterile.
- Do not use if device is damaged or expired.
- Not intended for permanent use. The Bluthane® Ureteral Stent must not remain indwelled for more than 180 days. If appropriate for the patient, the stent may be replaced with a new stent.
- Ensure compatibility of the guidewire before use. The compatible guidewire information is mentioned the device label.
- Product use should not be carried out by unauthorized or untrained persons. Product should be used by trained physician with support from qualified technicians using aseptic technique.
- On completion of the procedure, dispose off the device and its components, if any, as per institutional and local regulatory guidelines for bio-hazardous medical waste disposal.
- Keep away from fluorescent light, electrical motors and generators.
- Avoid contact with sharp and abrasive objects.
- Use only accessories supplied with package.
- Bluthane® Ureteral Stent can be removed by gently withdrawing using endoscopic forceps. Do not force if resistance is encountered. Device removal should be carried out by trained physician with support from qualified technicians.

PRODUCT RECOMMENDATIONS

If no wire guide is provided in the device package, the following is recommended.

Stent size (in FR)	Accepted wire guide size (in Inch)
2.8, 3, 3.5 & 3.8	0.018
4	0.025
4.5, 4.7 & 4.8	0.032
5, 5.5, 6 & 6.5	0.035
7 & 8	0.038

SUPPLY AND STORAGE

The device is packed with primary pack made of polyethylene and the secondary package is made of polyethylene and Tyvek pouch. Supplied sterilized by Ethylene Oxide gas in peel open package.

Device shall be stored at temperature range from 20°C to 40°C and away from direct sunlight and other adverse environmental conditions.

Bluthane® Ureteral Stents and Stent Sets

INSTRUCTIONS FOR USE

NOTE: Prior to use, If stent is hydrophilic coated, immerse stent in sterile water or isotonic saline to allow the hydrophilic coat activated and become lubricious.

This will allow the easy placement of stent under normal conditions.

BOTH END OPEN-TYPE:

1. Determine the appropriate stent length as per the patient requirement. This is generally calculated from the baseline pyelogram. A proper stent length enhances patient comfort and efficiency.
2. Insert soft tip of the compatible guidewire into the endoscope, access ureteral orifice and pass the guidewire up to the renal pelvis.
3. Pass tapered & un-marked tip of stent over guidewire and through endoscope while maintaining guidewire position.
4. Insert the pusher to the proximal end of the guide wire and advance stent up to the renal pelvis using the pusher under direct vision.
5. Observe for the distal end of the stent to appear at the ureterovesical junction. Once correct placement within the renal pelvis is confirmed by fluoroscopic guidance.
6. Remove the guidewire and then the pusher carefully. The stent pigtail loop will form once the guidewire is removed
7. If necessary, fine adjustments can be made using endoscopic forceps or with string if provided along with stent.
8. Detach the monofilament string from the stent and remove from patient body after stent adjustment.

ONE END CLOSED-TYPE:

1. Determine the proper stent length for the patient. This is generally calculated from the baseline pyelogram. A proper stent length enhances patient comfort and efficiency.
2. Insert the compatible guidewire with stiff tip into the stent and stretch the stent's coils at both ends.
3. Insert and connect the pusher to the stent from the guidewires proximal end and fix it to meet the cross section of the stent. Use the clamp to secure the pusher with the guidewire if necessary
4. Advance the stent, guidewire and pusher assembly until the tip reaches the renal pelvis under direct vision.
5. Observe for the distal end of the stent to appear at the ureterovesical junction. Confirm correct placement within the renal pelvis by fluoroscopic guidance.
6. Remove the guidewire and then the pusher carefully.
7. If necessary, fine adjustments can be made using endoscopic forceps or with string if provided along with stent.
8. Detach the string from the stent and remove from patient body after stent adjustment.

Summary of Safety and Clinical Performance (SSCP)

Summary of Safety and Clinical Performance for the device is available on the website of EUDAMED:<https://ec.europa.eu/tools/eudamed>

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
Non-pyrogenic

Apirógeno
Apirogeno

Single sterile barrier system with protective package inside

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Blueneem, the Blueneem Logo are Registered and Trademarks of Blue Neem Medical Devices Private Limited Company.

© 2024 All rights reserved

BLUE NEEM MEDICAL DEVICES PRIVATE LIMITED

Plot Nos 270 & 271, Road No 5,
Harohalli Industrial Area, II Phase,
Kanakapura Taluk, Ramanagara -
562112, India.

Phone: +91 9739972855 / 54

E-Mail: contact@blueneem.com

Website: www.blueneem.com

**Authorized Representative in the
European Community:**

Arazy Group (Ireland) Ltd

19 Baggot Street Lower,
Dublin 2, D02 X658, Ireland

E-mail: Ireland@arazygroup.com