# **Fascial Dilators and Dilator Sets**

## **Instructions for Use**

**EN**

**Fascial Dilators and Dilator Sets**

Instructions for use

**ES**

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# **Fascial Dilators and Dilator Sets**

## **DEVICE DESCRIPTION**

Fascial Dilators are tube like device made of medical grade plastic material, with a long slender hollow body and a tapered tip at the distal end. They are designed for stretching or enlarging the fascial tissue covering a cavity, tract or opening during a procedure. The Dilators are radiopaque for imaging guidance. Device is available in individual dilators (contents may varies and dilator sets (contents may varies based on requirement )

## **INTENDED PURPOSE & DURATION**

Fascial Dilators are used for fascial dilation, intended for stretching or enlarging the fascial tissue covering a cavity, tract or opening prior to a diagnostic or interventional procedure. up to a period of 60 minutes.

## **INDICATIONS FOR USE & CONTRAINDICATIONS**

None known.

## **POTENTIAL COMPLICATIONS**

Not limited to Urinary tract infection & Hemorrhage.

## **WARNINGS AND PRECAUTIONS**

- Use with caution when advancing the dilators to avoid tissues and orifice damages. Use only compatible wire guide.
- Do not use if product is damaged or expired, reuse and resterlization of this device leads to loss in structural integrity and performance.
- This device is intended for use by trained physicians and this device should not be used without comprehensive knowledge of the indications, techniques and risks of the procedure.
- On completion of the procedure, dispose off the device as per local regulatory guidelines for bio-hazardous medical waste disposal.

## **PROCEDURE**

1. Prior to use, if the product supplied with hydrophilic coating, immerse the device in sterile water or isotonic solution which makes it lubricious under normal conditions of use.
2. Perform initial access with puncture needle for the purpose of sequential dilation using standard procedure.
3. Introduce the flexible tip of the guidewire into needle which has previously been positioned into the renal collecting system under fluoroscopic control.
4. Advance the guidewire to the desired location. Once the guidewire is in position, remove the puncture needle from the wire guide prior to introduction of the dilators over the guidewire.
5. Progressing from the smallest to the largest appropriate size, pass the radiopaque dilators over the guidewire while maintaining the guidewire's position.

# **Fascial Dilators and Dilator Sets**

## **SUPPLY AND STORAGE**

The device is packed with primary pack made of polyethylene and the secondary package is made of polyethylene and Tyvek pouch. Supplied sterilized by Ethylene Oxide gas in peel open package.

Device shall be stored at temperature range from 20<sup>0</sup>C to 40<sup>0</sup>C and away from direct sunlight and other adverse environmental conditions.

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
**Non-pyrogenic**

Apirógeno  
Apirogeno

**Single sterile barrier system with protective package inside**

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Blueneem, the Blueneem Logo are Registered and Trademarks of Blue Neem Medical Devices Private Limited Company.

© 2024 All rights reserved

## **BLUE NEEM MEDICAL DEVICES PRIVATE LIMITED**

Plot Nos 270 & 271, Road No 5,  
Harohalli Industrial Area, II Phase,  
Kanakapura Taluk, Ramanagara -  
562112, India.

Phone: +91 9739972855 / 54

E-Mail: [contact@blueneem.com](mailto:contact@blueneem.com)

Website: [www.blueneem.com](http://www.blueneem.com)

---

**Authorized Representative in the  
European Community:**

**Arazy Group (Ireland) Ltd**

19 Baggot Street Lower,  
Dublin 2, D02 X658, Ireland

E-mail: [Ireland@arazygroup.com](mailto:Ireland@arazygroup.com)