# **Malecot Catheters**

## **Instructions for Use**

**EN**      **Malecot Catheters**  
Instructions for use

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# Malecot Catheters

## DEVICE DESCRIPTION

Malecot Catheter is flexible tube made of polyurethane, intended to be inserted into but not limited into the renal pelvis percutaneously and to gain access into the upper urinary tract and bladder percutaneously. The tip of the device has a wings shaped configuration which enhances retention. Device packaging is provided with plastic/metal retention straightener, urine bag connector. For product specifications refer to the label.

Malecot Catheters package includes(package components may vary):

- Radiopaque catheter (with or without Hydrophilic coating)
- Plastic/metal cannulated straightener
- Urine bag connector

## INTENDED PURPOSE & DURATION

Malecot catheters are intended to be inserted percutaneously into the renal pelvis and to gain access to the upper urinary tract to perform external drainage of urine from the renal pelvis for a short-term period of up to 30 days.

## INDICATIONS FOR USE & CONTRAINDICATIONS

Indicated when diagnosed to place drainage catheter percutaneously in kidney for renal drainage.

No known contraindications.

Device should not be used in the conditions of Bleeding Diathesis and uncontrolled hypertension, Anticoagulant use.

## POTENTIAL COMPLICATIONS

Potential complications include but not limited to Urinary tract infection, perforation of Kidney.

## WARNINGS AND PRECAUTIONS

- If a catheter has become malpositioned or if drainage ceases, the catheter should be promptly exchanged or removed.
- Do not use if product is damaged or expired and do not reuse the device.
- Use only accessories supplied along with package and do not resterilize.
- Manipulation of the product needs ultrasound or fluoroscopic guidance.
- If the device is hydrophilic coated, activate the coat by using sterile water or saline or other isotonic solution to make the catheter insertion easier.
- Dispose off the device and its components, if any, as per institutional and local regulatory guidelines for bio-hazardous medical waste disposal
- This device is intended for use by physicians trained and experienced in diagnostic and interventional techniques, standard techniques for placement of percutaneous drainage catheters should be employed.
- Periodic evaluation of the catheter is advised; catheter should not remain indwelled for more than 30 days after deployment.

# Malecot Catheters

## PRODUCT RECOMMENDATIONS

For catheter loading Seldinger technique is employed. Recommended wire guide size is mentioned as follows:

Catheter size ( in FR )	Accepted wire guide size ( in Inch )
8	0.035
10, 12, 14, 16, 18, 20 & 22	0.038

## INSTRUCTIONS FOR USE

### PROCEDURE: CATHETER PLACEMENT

1. Perform procedure using standard procedures of puncture, compatible guidewire placement and dilatation.
2. As required, Use the Plastic or Metal Straightener provided to close the Malecot wings and twist the hub to lock. Then pass the wings end of the catheter over the external end of the guidewire.
3. Gradually advance the wings end well into the collecting system.
4. Confirm the position using fluoroscopic guidance.
5. Unlock the Straightener with a gentle twist and remove the Straightener from the catheter to open the Malecot wings.
6. Using the Urine bag connector, connect the catheter to the urine collection bag.

### PROCEDURE: CATHETER REMOVAL

1. Detach the catheter from catheter fixing disc or skin.
2. Gently pull the catheter carefully to remove.
3. If necessary, place a guidewire, into the kidney/collecting system to maintain access.

## SUPPLY AND STORAGE

The device is packed with primary pack made of polyethylene and the secondary package is made of polyethylene and Tyvek pouch. Supplied sterilized by Ethylene Oxide gas in peel open package.

Device shall be stored at temperature range from 20<sup>0</sup>C to 40<sup>0</sup>C and away from direct sunlight and other adverse environmental conditions.

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
**Non-pyrogenic**

Apirógeno  
Apirogeno

**Single sterile barrier system with protective package inside**

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Blueneem, the Blueneem Logo are Registered and Trademarks of Blue Neem Medical Devices Private Limited Company.

© 2024 All rights reserved

## **BLUE NEEM MEDICAL DEVICES PRIVATE LIMITED**

Plot Nos 270 & 271, Road No 5,  
Harohalli Industrial Area, II Phase,  
Kanakapura Taluk, Ramanagara -  
562112, India.

Phone: +91 9739972855 / 54

E-Mail: [contact@blueneem.com](mailto:contact@blueneem.com)

Website: [www.blueneem.com](http://www.blueneem.com)

---

**Authorized Representative in the  
European Community:**

**Arazy Group (Ireland) Ltd**

19 Baggot Street Lower,  
Dublin 2, D02 X658, Ireland

E-mail: [Ireland@arazygroup.com](mailto:Ireland@arazygroup.com)