# **Mono J Catheter and Catheter sets**

## **Instructions for Use**

**EN**      **Mono J Catheters and Catheter Sets**  
Instructions for use

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# **Mono J Catheter and Catheter sets**

## **DEVICE DESCRIPTION**

Mono J Catheter is a flexible tube-like device made from polyurethane for introduction into urinary tract through an endoscope. Device is available in variety of French sizes, pigtail distal tip configurations for retention. Choice of device is based on clinical situation. Refer to label for technical specifications.

Device package includes (contents may varies):

- Radiopaque pigtail catheter
- Urine bag connector
- PTFE Wire guide ( available only in catheter sets)

## **INTENDED PURPOSE & DURATION**

Used for temporary drainage and irrigation of urine from the kidney for a period of up to 48 hours.

## **INDICATIONS FOR USE & CONTRAINDICATIONS**

Contraindicated for use in patients unsuitable for catheterization.

## **POTENTIAL COMPLICATIONS**

Potential complications Include but not limited to hemorrhage, edema, perforation of the kidney, ureter avulsion renal pelvis, ureter and/or bladder, urinary tract infection.

## **WARNINGS AND PRECAUTIONS**

- Manipulation of device requires Ultrasound / Fluoroscopy / other image guiding techniques.
- Do not use if product is damaged or expired, reuse and resterilization of this device leads to potential complications and loss in structural integrity
- This device is intended for use by physicians trained and experienced in diagnostic and endoscopic techniques, standard techniques for placement of endoscopic devices should be employed.
- On completion of the procedure, dispose off the device and its components, if any, as per institutional and local regulatory guidelines for bio-hazardous medical waste disposal.
- Not intended for permanent use.
- Use only compatible endoscopic devices.
- Avoid bending or kinking of catheter prior to placement. Doing so could damage the integrity of the catheter and result in patient injury.

# **Ureteral Catheter**

## **INSTRUCTIONS FOR USE**

### **PROCEDURE: FOR CLOSED END.**

1. Insert the compatible guidewire with stiff tip into the straight end of the catheter and stretch the catheter pigtail coil at the other end.
2. Advance the catheter and guidewire until the tip reaches the renal pelvis.
3. Once correct placement within the renal pelvis is confirmed by fluoroscopic guidance, remove the guidewire carefully.

### **PROCEDURE: FOR OPEN END**

1. Insert soft tip of the compatible guidewire into endoscope, access ureteral orifice and pass the guidewire up to the renal pelvis
2. Pass the catheter over guidewire, loop first, through the endoscope and while maintaining guidewire position.
3. Push the catheter until it reaches renal pelvis under fluoroscopic guidance.
4. Once correct placement within the renal pelvis is confirmed by fluoroscopic guidance, remove the guidewire carefully. Remove the endoscope while maintaining the position of the catheter.

## **PRODUCT RECOMMENDATIONS**

Recommended wire guide for the usage of this product is as follows:

4.7 Fr : 0.032 inch wire guide ; 05 Fr : 0.035 inch wire guide ; 06, 07 & 08 Fr : 0.038 inch wire guide.

## **SUPPLY AND STORAGE**

The device is packed with primary pack made of polyethylene and the secondary package is made of polyethylene and Tyvek pouch. Supplied sterilized by Ethylene Oxide gas in peel open package.

Device shall be stored at temperature range from 20<sup>0</sup>C to 40<sup>0</sup>C and away from direct sunlight and other adverse environmental conditions.

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
**Non-pyrogenic**

Apirógeno  
Apirogeno

**Single sterile barrier system with protective package inside**

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Blueneem, the Blueneem Logo are Registered and Trademarks of Blue Neem Medical Devices Private Limited Company.

© 2024 All rights reserved

## **BLUE NEEM MEDICAL DEVICES PRIVATE LIMITED**

Plot Nos 270 & 271, Road No 5,  
Harohalli Industrial Area, II Phase,  
Kanakapura Taluk, Ramanagara -  
562112, India.

Phone: +91 9739972855 / 54

E-Mail: [contact@blueneem.com](mailto:contact@blueneem.com)

Website: [www.blueneem.com](http://www.blueneem.com)

---

**Authorized Representative in the  
European Community:**

**Arazy Group (Ireland) Ltd**

19 Baggot Street Lower,  
Dublin 2, D02 X658, Ireland

E-mail: [Ireland@arazygroup.com](mailto:Ireland@arazygroup.com)