---

---

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# **Percutaneous Nephrostomy Drainage Set - Malecot**

## **DEVICE DESCRIPTION**

Percutaneous Nephrostomy Drainage Catheter - malecot is flexible tube made of polyurethane, intended to be inserted into the renal pelvis percutaneously and to gain access into the upper urinary tract. The tip of the device has wings shaped configuration which enhances retention. For product specifications refer to the label.

Percutaneous Nephrostomy Drainage Set – Malecot package includes(set components may vary):

- Radiopaque malecot catheter (with or without Hydrophilic coating)
- Plastic/metal cannulated loop straightener
- PTFE Guidewire
- Initial puncture needle
- Catheter fixing disc
- Scalpel blade
- Fascial dilators
- Urine bag connector

## **INTENDED PURPOSE & DURATION**

Percutaneous Nephrostomy Drainage Set “Malecot” is intended for percutaneous placement of pigtail catheter into the renal pelvis, to gain access to the upper urinary tract to perform external drainage of urine from the renal pelvis for a short-term period of up to 30 days.

## **INDICATIONS FOR USE & CONTRAINDICATIONS**

No known contraindications.

Device should not be used in the conditions of Bleeding Diathesis and uncontrolled hypertension, Anticoagulant use.

## **POTENTIAL COMPLICATIONS**

Malposition, Infection, Catheter dislodgement and pain.

## **WARNINGS AND PRECAUTIONS**

- If a catheter has become malpositioned or if drainage ceases, the catheter should be promptly exchanged or removed.
- Do not use if product is damaged or expired and do not reuse the device.
- Use only accessories supplied along with package to avoid compatibility issue.
- Manipulation of the product needs ultrasound or fluoroscopic guidance.
- If the device is hydrophilic coated, activate the coat by using sterile water or saline or other isotonic solution to make the catheter insertion easier.

# **Percutaneous Nephrostomy Drainage Set - Malecot**

- Retain the plastic/metal cannulated straightener to straighten the catheter loop to remove the catheter.
- Dispose off the device and its components, if any, as per institutional and local regulatory guidelines for bio-hazardous medical waste disposal
- This device is intended for use by physicians trained and experienced in diagnostic and interventional techniques, standard techniques for placement of percutaneous drainage catheters should be employed.
- Periodic evaluation of the catheter is advised; catheter should not remain indwelled for more than 30 days after deployment.
- Do not forcefully advance any components during removal or placement of the catheter. Carefully remove the components if any resistance is occurred.
- Reuse and Resterlization of the device is strictly prohibited.

## **INSTRUCTIONS FOR USE**

NOTE: Prior to use, If catheter is hydrophilic coated, immerse stent in sterile water or isotonic saline to allow the hydrophilic coat activated and become lubricious. This will allow the easy placement of stent under normal conditions.

- By retrograde pyelogram, ultrasound or CT scan, localize the kidney/collection system to be drained.
- After pyelogram, Identify and mark the skin site overlaying the kidney/collection system. Anesthetize the marked skin site.
- With the help of scalpel, make a small incision.
- Under fluoroscopic guidance, introduce the IP needle provided along with set, at right angles to the skin at the marked site. After achieving placement, remove the inner trocar of the needle. Confirm placement via free flow of urine.
- On confirmation, pass the guidewire through the outer sheath of the IP Needle.
- On reaching the desired position, remove the outer sheath of the needle while maintaining guidewire position.
- To facilitate passage of the catheter, dilate the musculofascial track using dilators in sequence from smallest to the largest.
- Use the Straightner provided to close the Malecot wings and twist the hub to lock. Then pass the wings end of the catheter over the external end of the guidewire. Gradually advance the wings end well into the collecting system. Confirm the position using fluoroscopic guidance.
- Unlock the Straightner with a gentle twist and remove the Straightner from the catheter to open the Malecot wings.
- Holding the shaft of the catheter securely in position remove the guidewire carefully.

# **Percutaneous Nephrostomy Drainage Set - Malecot**

- Tape or suture the catheter fixing disc to the skin.
- Secure the catheter to the catheter fixing disc.
- Using the urine bag connector, connect the catheter to the urine bag.

## **PROCEDURE: CATHETER REMOVAL FOR STANDARD TYPE**

- Detach the catheter from catheter fixing disc.
- Gently pull the catheter carefully to remove.
- If necessary, place a guidewire, into the kidney/collecting system to maintain access.

## **SUPPLY AND STORAGE**

The device is packed with primary pack made of polyethylene and the secondary package is made of polyethylene and Tyvek pouch. Supplied sterilized by Ethylene Oxide gas in peel open package.

Device shall be stored at temperature range from 20°C to 40°C and away from direct sunlight and other adverse environmental conditions.

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
**Non-pyrogenic**

Apirógeno  
Apirogeno

**Single sterile barrier system with protective package inside**

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Blueneem, the Blueneem Logo are Registered and Trademarks of Blue Neem Medical Devices Private Limited Company.

© 2024 All rights reserved

## **BLUE NEEM MEDICAL DEVICES PRIVATE LIMITED**

Plot Nos 270 & 271, Road No 5,  
Harohalli Industrial Area, II Phase,  
Kanakapura Taluk, Ramanagara -  
562112, India.

Phone: +91 9739972855 / 54

E-Mail: [contact@blueneem.com](mailto:contact@blueneem.com)

Website: [www.blueneem.com](http://www.blueneem.com)

---

**Authorized Representative in the  
European Community:**

**Arazy Group (Ireland) Ltd**

19 Baggot Street Lower,  
Dublin 2, D02 X658, Ireland

E-mail: [Ireland@arazygroup.com](mailto:Ireland@arazygroup.com)