Re-Entry Malecot Catheters

Instructions for Use

EN **Re-Entry Malecot Catheters**
Instructions for use

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Re-Entry Malecot Catheters

DEVICE DESCRIPTION

Re-Entry Malecot Catheter is flexible tube made of polyurethane, intended to be inserted into the renal pelvis percutaneously with self-retaining Malecot winged tip to perform external drainage of urine from kidney and long tapered end inserted into ureter to perform internal drainage to bladder. Device packaging is provided with two plastic/metal retention straighteners, urine bag connector. For product specifications refer to the label.

Malecot Catheters package includes(package components may vary):

- Radiopaque catheter (with or without Hydrophilic coating)
- Plastic/metal cannulated straighteners
- Urine bag connector

INTENDED PURPOSE & DURATION

Re-Entry Malecot Drainage Catheter is intended to be inserted percutaneously into the renal pelvis and to gain access to the upper urinary tract & the ureter to perform external drainage of urine from the renal pelvis and internal drainage into the bladder.

INDICATIONS FOR USE & CONTRAINDICATIONS

Indicated when diagnosed to place drainage catheter percutaneously in kidney for renal drainage. No known contraindications.

Device should not be used in the conditions of Bleeding Diathesis and uncontrolled hypertension, Anticoagulant use.

POTENTIAL COMPLICATIONS

Potential complications include catheter occlusion and/or dislodgement, hemorrhage or hematoma, infection and pain.

WARNINGS AND PRECAUTIONS

- If a catheter has become malpositioned or if drainage ceases, the catheter should be promptly exchanged or removed.
- Do not use if product is damaged or expired and do not reuse the device.
- Use only accessories supplied along with package and do not resterilize.
- Manipulation of the product needs ultrasound or fluoroscopic guidance.
- If the device is hydrophilic coated, activate the coat by using sterile water or saline or other isotonic solution to make the catheter insertion easier.
- Dispose off the device and its components, if any, as per institutional and local regulatory guidelines for bio-hazardous medical waste disposal.
- This device is intended for use by physicians trained and experienced in diagnostic and interventional techniques, standard techniques for placement of percutaneous drainage catheters should be employed.
- Periodic evaluation of the catheter is advised; catheter should not remain indwelled for more than 30 days after deployment.

Malecot Catheters

PRODUCT RECOMMENDATIONS

For catheter loading use 0.035 inch wire guide.

INSTRUCTIONS FOR USE

PROCEDURE: CATHETER PLACEMENT

1. Perform the standard techniques for establishing a nephrostomy tract utilizing fluoroscopic control by leaving a 150 cm, 0.035 Inch guidewire in the collecting system.
2. Confirm via fluoroscopy that the guidewire wire travels through all the length of the ureter and is well into the bladder.
3. Dilate the musculofascial tract by progressing dilators sequentially from smallest size to the largest size to facilitate the passage of the catheter.
4. Insert the Straightener provided along with the catheter and gently twist the hub to lock and close the wings and straighten the catheter.
5. Pass the thin end of the catheter over the external end of the guidewire and gradually advance the catheter wings well into the collecting system and place the thin end of the catheter into the ureter.
6. Confirm the position fluoroscopically then release the Straightener from the catheter to open Malecot flower in renal pelvis. Carefully remove the Straightener.
7. Remove the guidewire holding the catheter in place and secure the catheter to the skin or to the fixing disc.
8. Connect the catheter to the drainage bag using the urine bag connector.

PROCEDURE: CATHETER REMOVAL

1. Detach the catheter from catheter fixing disc.
2. Gently pull the catheter carefully to remove.
3. If necessary, place a guidewire, into the kidney/collecting system to maintain access.

SUPPLY AND STORAGE

The device is packed with primary pack made of polyethylene and the secondary package is made of polyethylene and Tyvek pouch. Supplied sterilized by Ethylene Oxide gas in peel open package.

Device shall be stored at temperature range from 20⁰C to 40⁰C and away from direct sunlight and other adverse environmental conditions.

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
Non-pyrogenic

Apirógeno
Apirogeno

Single sterile barrier system with protective package inside

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Blueneem, the Blueneem Logo are Registered and Trademarks of Blue Neem Medical Devices Private Limited Company.

© 2024 All rights reserved

BLUE NEEM MEDICAL DEVICES PRIVATE LIMITED

Plot Nos 270 & 271, Road No 5,
Harohalli Industrial Area, II Phase,
Kanakapura Taluk, Ramanagara -
562112, India.

Phone: +91 9739972855 / 54

E-Mail: contact@blueneem.com

Website: www.blueneem.com

**Authorized Representative in the
European Community:**

Arazy Group (Ireland) Ltd

19 Baggot Street Lower,
Dublin 2, D02 X658, Ireland

E-mail: Ireland@arazygroup.com