Renal Dilators

Instructions for Use

EN

Renal Dilators

Instructions for use

ES

Dilatadores renales

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IT

Dilatatori renali

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Renal Dilators

DEVICE DESCRIPTION

Renal Dilators are tube like device made of medical grade plastic material, with a long slender hollow body and a tapered tip at the distal end. They are designed for stretching or enlarging the fascial tissue covering a cavity, tract or opening during a procedure. The Dilators are radiopaque for fluoroscopic guidance.

INTENDED PURPOSE & DURATION

Used for dilation of a percutaneous tract and setup of a working channel into the kidney for removal of stones up to a period of 60 minutes.

INDICATIONS FOR USE & CONTRAINDICATIONS

None known

POTENTIAL COMPLICATIONS

Not limited to Urinary tract infection & Hemorrhage.

WARNINGS AND PRECAUTIONS

- Use with caution when advancing the sheath and dilators to avoid tissues and orifice damages. Use only compatible advancing sheath.
- Do not use if product is damaged or expired, reuse and reesterlization of this device leads to loss in structural integrity and performance.
- This device is intended for use by trained physicians and this device should not be used without comprehensive knowledge of the indications, techniques and risks of the procedure.
- On completion of the procedure, dispose off the device as per local regulatory guidelines for bio-hazardous medical waste disposal.

PROCEDURE

1. Perform initial access with puncture needle, place a guidewire and perform sequential dilation as appropriate using standard procedure.
2. Now introduce the dilator of required size. Advance to the desired location using imaging guidance.
3. Introduce the Amplatz sheath over the dilator, use the amplatz sheath of appropriate size or use the same amplatz sheath if supplied along with the dilator.
4. Remove the dilator, leaving the amplatz sheath in position.
5. This can now be used as a working channel to perform the procedure.

SUPPLY AND STORAGE

The device is packed with primary pack made of polyethylene and the secondary package is made of polyethylene and Tyvek pouch. Supplied sterilized by Ethylene Oxide gas in peel open package.

Device shall be stored at temperature range from 20^oC to 40^oC and away from direct sunlight and other adverse environmental conditions.

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
Non-pyrogenic

Apirógeno
Apirogeno

Single sterile barrier system with protective package inside

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Blueneem, the Blueneem Logo are Registered and Trademarks of Blue Neem Medical Devices Private Limited Company.

© 2024 All rights reserved

BLUE NEEM MEDICAL DEVICES PRIVATE LIMITED

Plot Nos 270 & 271, Road No 5,
Harohalli Industrial Area, II Phase,
Kanakapura Taluk, Ramanagara -
562112, India.

Phone: +91 9739972855 / 54

E-Mail: contact@blueneem.com

Website: www.blueneem.com

**Authorized Representative in the
European Community:**

Arazy Group (Ireland) Ltd

19 Baggot Street Lower,
Dublin 2, D02 X658, Ireland

E-mail: Ireland@arazygroup.com