Single Step Ureteral Dilators

Instructions for Use

EN

Single Step Ureteral Dilators

Instructions for use

ES

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Single Step Ureteral Dilators

DEVICE DESCRIPTION

Single Step Ureteral Dilators are long tube like device made of medical grade plastic material, with a long tapered tip at the distal end. The Dilators are radiopaque for imaging guidance.

INTENDED PURPOSE & DURATION

Single Step Ureteral Dilator is used for one step dilation of the ureter through natural orifice prior to a diagnostic or interventional procedure up to a period of 04 hours.

INDICATIONS FOR USE & CONTRAINDICATIONS

None known.

POTENTIAL COMPLICATIONS

Edema, Irritation, inflammation, Ureteral avulsion.

WARNINGS AND PRECAUTIONS

- Use with caution when advancing the dilators to avoid tissues and orifice damages. Use only compatible wire guide.
- Do not use if product is damaged or expired, reuse and resterilization of this device leads to loss in structural integrity and performance.
- This device is intended for use by trained physicians and this device should not be used without comprehensive knowledge of the indications, techniques and risks of the procedure.
- On completion of the procedure, dispose off the device as per local regulatory guidelines for bio-hazardous medical waste disposal.

PROCEDURE

1. Activate the hydrophilic coating by using sterile aqueous media before use, if dilator is supplied with hydrophilic coating.
2. Under fluoroscopic guidance, pass the guidewire to the desired distance into the ureter. Insert the dilator on to the guidewire
3. Progressing from the tip to the largest appropriate size tapered on the dilator; pass the radiopaque dilator over the guidewire while maintaining the Guidewire's position to the desired distance.

SUPPLY AND STORAGE

The device is packed with primary pack made of polyethylene and the secondary package is made of polyethylene and Tyvek pouch. Supplied sterilized by Ethylene Oxide gas in peel open package.

Device shall be stored at temperature range from 20°C to 40°C and away from direct sunlight and other adverse environmental conditions.

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
Non-pyrogenic

Apirógeno
Apirogeno

Single sterile barrier system with protective package inside

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Blueneem, the Blueneem Logo are Registered and Trademarks of Blue Neem Medical Devices Private Limited Company.

© 2024 All rights reserved

BLUE NEEM MEDICAL DEVICES PRIVATE LIMITED

Plot Nos 270 & 271, Road No 5,
Harohalli Industrial Area, II Phase,
Kanakapura Taluk, Ramanagara -
562112, India.

Phone: +91 9739972855 / 54

E-Mail: contact@blueneem.com

Website: www.blueneem.com

**Authorized Representative in the
European Community:**

Arazy Group (Ireland) Ltd

19 Baggot Street Lower,
Dublin 2, D02 X658, Ireland

E-mail: Ireland@arazygroup.com