Ureteral Catheter

Instructions for Use

EN **Ureteral Catheters**
Instructions for use

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Ureteral Catheter

DEVICE DESCRIPTION

Ureteral Catheter is a flexible tube-like device made from PEBA for introduction into urinary tract through an endoscope. Device is available in variety of French sizes, lengths and distal tip configurations. Choice of device is based on clinical situation. Refer to label for technical specifications.

INTENDED PURPOSE & DURATION

Ureteral Catheters are used for catheterization of the urinary tract during interoperative procedures, includes the following operations:

- Injection of contrast media such as to perform retrograde pyelogram.
- Navigation of a tortuous ureter.
- Access, drainage of fluids and advancement or exchange of ancillary devices including guidewires (for open end catheters only) for a period of up to 06 hours.

INDICATIONS FOR USE & CONTRAINDICATIONS

Contraindicated for use in patients unsuitable for catheterization.

POTENTIAL COMPLICATIONS

Potential complications Include but not limited to hemorrhage, edema, perforation of the kidney, ureter avulsion renal pelvis, ureter and/or bladder, urinary tract infection.

WARNINGS AND PRECAUTIONS

- Manipulation of device requires Ultrasound / Fluoroscopy / Other image guiding techniques.
- Do not use if product is damaged or expired, reuse and resterilization of this device leads to potential complications and loss in structural integrity
- This device is intended for use by physicians trained and experienced in diagnostic and endoscopic techniques, standard techniques for placement of endoscopic devices should be employed.
- On completion of the procedure, dispose off the device and its components, if any, as per institutional and local regulatory guidelines for bio-hazardous medical waste disposal.
- Not intended for permanent use.
- Use only compatible endoscopic devices.
- Avoid bending or kinking of catheter prior to placement. Doing so could damage the integrity of the catheter and result in patient injury.

Ureteral Catheter

INSTRUCTIONS FOR USE

PROCEDURE: FOR CLOSED END.

Supplied with Stylet (Placed without the guidewire)

1. Place the endoscope into bladder and view the ureteral orifice, Choose the appropriate catheter type based on patient anatomy.
2. The catheter adapter is secured to the ureteral catheter.
3. Ensure that the stylet is inserted into the catheter.
4. Under fluoroscopic vision the catheter tip is positioned within the ureteral orifice through the endoscope.
5. Remove the stylet carefully.
6. Contrast media can then be injected through the ureteral catheter into the ureter.

PROCEDURE: FOR OPEN END

(Placed with the guidewire)

1. Place the endoscope into the bladder and view the ureteral orifice. Choose the appropriate catheter type based on patient anatomy.
2. Under fluoroscopic guidance, place the guidewire within the ureteral orifice.
3. The ureteral catheter is inserted over the guidewire and advanced into the ureter.
4. Remove the guidewire carefully.
5. Secure the adapter to the catheter.
6. Contrast media can then be injected through the ureteral catheter into the ureter.

PRODUCT RECOMMENDATIONS

Recommended wire guide for the usage of this product is as follows:

03 Fr : 0.018 inch wire guide ; 04 Fr : 0.025 inch wire guide ; 05 Fr : 0.035 inch wire guide ; 06, 07 & 08 Fr : 0.038 inch wire guide.

SUPPLY AND STORAGE

The device is packed with primary pack made of polyethylene and the secondary package is made of polyethylene and Tyvek pouch. Supplied sterilized by Ethylene Oxide gas in peel open package.

Device shall be stored at temperature range from 20⁰C to 40⁰C and away from direct sunlight and other adverse environmental conditions.

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
Non-pyrogenic

Apirógeno
Apirogeno

Single sterile barrier system with protective package inside

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Blueneem, the Blueneem Logo are Registered and Trademarks of Blue Neem Medical Devices Private Limited Company.

© 2024 All rights reserved

BLUE NEEM MEDICAL DEVICES PRIVATE LIMITED

Plot Nos 270 & 271, Road No 5,
Harohalli Industrial Area, II Phase,
Kanakapura Taluk, Ramanagara -
562112, India.

Phone: +91 9739972855 / 54

E-Mail: contact@blueneem.com

Website: www.blueneem.com

**Authorized Representative in the
European Community:**

Arazy Group (Ireland) Ltd

19 Baggot Street Lower,
Dublin 2, D02 X658, Ireland

E-mail: Ireland@arazygroup.com