Hydro Twister® Guidewires

Instructions for Use

EN

Hydro Twister® Guidewires

Instructions for use

ES

Guías Hydro Twister®

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IT

Fili guida Hydro Twister®

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Hydro Twister® Guidewires

DEVICE DESCRIPTION

Hydro Twister® guidewires are made of a Nitinol core mandril inside and outer Polymer sleeve. Hydro Twister® guidewires are thin wires with flexible hydrophilic tips available in different diameters, lengths and straight flexible tip configurations that can be inserted into a tract to facilitate the placement of a medical device or instrument such as an endoscope, catheter, stent or dilator.

The surface of the wire is made in such a way to facilitate coaxial use with other devices. Please refer to the product label for product specifications.

Hydro Twister® Guidewires are available in following tip configuration:

- Straight tip

INTENDED PURPOSE & DURATION

Guidewire is used to gain access to desired anatomical location and to assist in placement of medical devices such as catheters, stents and dilators during diagnostic or interventional procedures for a short-term use of less than 04 hours.

INDICATIONS FOR USE & CONTRAINDICATIONS

None Known

POTENTIAL COMPLICATIONS

Include but not limited to Tissue perforation, Infection, Tissue trauma & Edema.

WARNINGS AND PRECAUTIONS

- Single use only.
- Do Not Reuse, reuse can cause severe infection. Intended use may fail on reuse and might cause kink of wire on reuse.
- Use on or before the shelf life printed on label.
- Do not resterilize & use only if packaging is intact.
- This device is not intended to be used for coronary arteries and other delicate tortuous anatomical regions.
- Avoid forceful angulation.

Hydro Twister® Guidewires

- Avoid manipulating or withdrawing the guidewire back through metal cannula or needle. A sharp edge may scrape or shear material from wire.
- Do not alter the tip configuration manually. It may damage the guidewire.
- Do not attempt to torque the guidewire.
- This device is intended to be used by trained physicians or trained technicians
- Use medical imaging when manipulating the guidewire. Do not advance or manipulate the guidewire without the visual evidence of the corresponding movement of the distal tip.
- When guidewire is used in association with another device, compatibility shall be checked before using in the procedure.
- Use only flexible end of the guidewire as initial inserted end. Inserting from rigid end may cause the damage to tissue or perforation.
- Inspect the wire for kinks, Sleeve defects or other damages prior to use.

INSTRUCTIONS FOR USE

- Open the device package and remove the guidewire using aseptic technique.
- Locate the flexible tip of the guidewire and ensure the guidewire compatibility with other device. Activate the hydrophilic tip by using isotonic solution then insert or advance the wire through the device or orifice and position the flexible tip of the guidewire in desired anatomical location.
- Use only medical imaging techniques for manipulating the guidewire.
- Perform the procedure and retrieve the guidewire gently.
- Discard the wire after the procedure.

SUPPLY AND STORAGE

The device is packed with dispenser and then packed in peel open Tyvek pouch made of polyethylene and Tyvek. Supplied sterilized by Ethylene Oxide gas in peel open package.

Device shall be stored at temperature range from 20⁰C to 40⁰C and away from direct sunlight and other adverse environmental conditions.

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Blueneem, the Blueneem Logo are Registered and Trademarks of Blue
© 2024 All rights reserved
Neem Medical Devices Private Limited Company.
Arazy Group (Ireland) Ltd