Dalela Triple J Ureteral Stent

Device Description

Dalela Triple J ureteral Stent are flexible tubular devices made from radiopaque polyurethane. One or both ends of the stent consist of curled pigtail loops. Both ends of the stent are open or one end of stent is closed based on the type chosen. Side holes are punched along the length of the stent including the pigtail loops of the stent. The stent is marked along the length to enable visualization, advancement and placement.

An additional 'J' like bend is created at bladder end, as a result of this bend the lower coil of stent gets directed towards the dome and away from the sensitive trigone, Material used - Advanced thermal sensitive polyurethane that softens the stent when indwelled. If required stent will be supplied with hydrophilic coating, this coating when activated will make stent outer surface slippery for easy insertion.

Dalela Triple J Ureteral Stent Package Includes

- Radiopaque curled loop stent
- Pusher/Stent Positioner with radiopaque tip
- Clamp

Dalela Triple J Ureteral Stent Set Package Includes

- Radiopaque curled loop stent
- Pusher/Stent Positioner with radiopaque tip
- Guidewire
- Clamp

Size (Fr)

3, 3.5, 3.7, 3.8, 4, 4.5, 4.7, 4.8, 5, 5.5, 6, 6.5, 7 & 8

Length (cm)

22, 24, 26 & 28

Variant

Both ends open

One end Closed

All the variants are available in with and without Hydrophilic coat

Intended Purpose

The Triple J Ureteral Stent is a non-bioabsorbable medical device intended to be implanted for treatment of obstructed ureter for drainage of urine from renal pelvis to urinary bladder and is intended for long-term use up to 180 Days.

Indications For Use & Contraindications

Indicated to be administered in an obstructed ureter (e.g., from stone, stricture, compression, fibrosis, trauma). Not to be implanted in conditions of extreme strictures caused by tumours.

Potential Complications

Include but not limited to Ureteral reflux, Urinary tract infection (Uro sepsis), Hematuria, Edema, Dysuria, Encrustation and stent fragmentation.

Warnings

- Do not reuse the device. Do not use if package is opened or damaged and if there is suspect that the device is not sterile.
- Do not use if device is expired.
- Use only accessories supplied with package.
- Device is not intended for permanent use. The Ureteral Stent must not remain indwelled for more than 180 days. If appropriate or required for the patient, the stent may be replaced with a new stent.
- Ensure guidewire compatibility before use. The compatible guidewire information is mentioned on the device label.
- Product use should not be carried out by unauthorized or untrained persons. Product should be used by trained physician with support from qualified technicians using aseptic technique. On completion of the procedure, dispose off the device and its components, if any, as per institutional and local regulatory guidelines for bio-hazardous medical waste disposal.

- Keep away from fluorescent light, electrical motors and generators.
- Store in a dark, cool and dry place.
- Avoid contact with sharp and abrasive objects.
- Pregnant patient must be carefully observed for possible stent encrustation because of calcium supplements.
- Improper handling or angulation of stent, positioner and guidewire can seriously weaken the stent. Acute bending or overstressing during placement may result in subsequent separation of the stent at the point of stress during indwelling period.
- Stents should be checked periodically for the signs of encrustation and proper function. Periodic checks at intervals deemed to be appropriate by the physician in consideration of the individual patient's condition and other patient specific factors. Besides, use of this device should be based upon risk-benefit factors as they apply to your patient.
- Migration of the Ureteral Stent is a possible complication, which could require medical intervention for removal.
- Do not force components during removal or replacement. Carefully remove the components if resistance is encountered.
- Care should be taken in selecting a stent of the right length depending on the patient anatomy.
- The use of the stent should be based on consideration of risk benefit analysis as it applies to your patient.
- ****The stent needs to be placed in reasonably full bladder. ****Lower coil must be placed towards the dome of the bladder.
- In most patients, however, as the bladder lumen is empty due to ongoing bladder drainage by a urethral catheter, and as the abdominal pressure in prone position pushes the dome of bladder downwards, the lower coil of stent may not be placed towards the dome. The stent at this point of time needs to be placed in reasonably full bladder. It is for this reason we recommend clamping of Foley urethral catheter in the beginning.

Product Recommendations

Size (Fr)	Accepted Guidewire (inch)
3, 3.5 & 3.8	0.018
4	0.025
4.5, 4.7 & 4.8	0.032
5, 5.5, 6 & 6.5	0.035
7 & 8	0.038

Instruction For Use

For **Stents supplied with Hydrophilic coated Guidewire:** Guidewire's hydrophilic property should be activated by injecting sterile water or isotonic saline into the dispenser tube and then retrieve the guidewire for better results.

Endoscopic Placement of the Device

Both Ends Open

Prior to use immerse the stent in sterile water or isotonic saline to activate the stent surface lubricious (for Hydrophilic coated stents only)

1. Determine the proper stent length for the patient. This is generally calculated from the baseline pyelogram. A proper stent length enhances patient comfort and efficiency.
2. Insert soft tip of the compatible guidewire into the endoscope, access ureteral orifice and pass the guidewire up to the renal pelvis.
3. Pass tapered & un marked tip of stent over guidewire and through endoscope while maintaining guidewire position.
4. Insert the pusher to the proximal end of the guidewire and advance stent up to the renal pelvis using the pusher under direct vision.

One End Closed

Prior to use immerse the stent in sterile water or isotonic saline to activate the stent surface lubricious (for Hydrophilic coated stents only)

1. Determine the proper stent length for the patient. This is generally calculated from the baseline pyelogram. A proper stent length enhances patient comfort and efficiency.
2. Insert the compatible guidewire with stiff tip into the stent and stretch the stent's coils at both ends.
3. Insert and connect the pusher to the stent from the guidewires proximal end (triple J end) and fix it to meet the cross section of the stent. Use the clamp to secure the pusher with the guidewire if necessary.
4. Advance the stent, guidewire and pusher assembly until the tip reaches the renal pelvis under direct vision.

5. Observe for the distal end of the stent to appear at the ureterovesical junction. Once correct placement within the renal pelvis is confirmed by fluoroscopic guidance.
6. Remove the guidewire and then the pusher carefully. The stent pigtail loop will form once the guidewire is removed.
7. If necessary, fine adjustments can be made using endoscopic forceps.

5. Observe for the distal end of the stent to appear at the ureterovesical junction. Confirm correct placement within the renal pelvis by fluoroscopic guidance.
6. Remove the guidewire and then the pusher carefully.
7. If necessary, fine adjustments can be made using endoscopic forceps.

Removal of the Device

The Triple J Ureteral Stent can be removed by gently withdrawing the stent using endoscopic forceps. Do not force if resistance is encountered. Carefully and gently withdraw the stent, Product removal should not be carried out by unauthorized or untrained persons. Product should be removed by trained physician with support from qualified technicians.

Target Population Group

Patients with ureteral injury for adults and pediatric patient groups and limitation of using appropriate size and length.

Report Any Serious Incident

If problems occur using this device or any serious incidents or adverse events, please contact our customer relations department or call to the contact provided in our manufacturing address and competent authority of the Member State in which the user and/or patient is established.

Supply and Storage

The device is packed in primary pack made of polyethylene and the secondary package made of polyethylene and Tyvek. Supplied sterilized by Ethylene Oxide gas in peel open package. Device shall be stored at temperature range from 20°C to 40°C and away from direct sunlight.

Summary on Safety and Clinical Performance

Summary of Safety and Clinical Performance of this device is available on the website
EUDAMED:<https://ec.europa.eu/tools/eudamed>

Symbol Glossory And Definations

Standard Reference:

ISO 15223-1 Medical devices — Symbols to be used with information to be supplied by the manufacturer

Symbol	Symbol Title
	Do not re-sterilize
	Do not use if package is damaged and consult instructions for use
	Consult instructions for use
	Temperature limit
	Keep away from sunlight
	Do not re-use
	Not made with natural rubber latex
	Sterilized using ethylene oxide
	Non-pyrogenic
	Keep dry
	Medical device
	Single sterile barriersystem with protective packaging inside
	Does not contain DEHP
	Unique device identifier
	Date of Manufacture
	Use-by date
	Batch code

	Catalogue number
	Authorized Representative in the European Community/ European Union
	Manufacturer
	Caution
	Prescription

- This device is to be sold by/on the directions of a physician. Sterile (ETO) if package is unopened or undamaged.
- Read all the instructions before use.
- On completion of the procedure, dispose off device and its components if any as per institutional and local regulatory guidelines for bio-hazardous medical waste dispose.

Authorized Representative in the European Community
Arazy Group (Ireland) Ltd
19 Baggot Street Lower,
Dublin 2, D02 X658, Ireland
E-mail: Ireland@arazygroup.com

Manufacturer
Blue Neem Medical Devices Pvt Ltd
Plot Nos 270 & 271, Road No 5, Harohalli Industrial Area
II Phase, Kanakapura Taluk, Ramanagara
Karnataka - 562112, India.
Phone: +91 80 29761336
E-Mail: contact@blueneem.com
Website: www.blueneem.com

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