Dual Lumen Ureteral Access Sheath

Device Description

Dual lumen Ureteral Access Sheath is a flexible tube-like instrument. It is designed to function as a barrier for the endoscope and other instruments from patient contact. It contains ports/ lumens for the insertion of instruments with wire guide backloading clip

Size (Fr)

9.5, 10, 11, 11.5 & 12

Length (cm)

28, 35, 40, 45 & 55

Intended Purpose

Dual Lumen Ureteral Access Sheath is used to establish a conduit during urological procedures facilitating the passage for Ureteroscope/endoscope and other surgical instruments into the urinary tract. The working channel protects the ureter during repeated instrument exchanges. The sheath also protects delicate instruments and smaller flexible scopes from damage.

Indications For Use & Contraindications

Indicated to be used in urinary obstruction due to urinary stones. Patients who are contraindicated for retrograde urological procedures; presence of tight strictures, which would limit use of the device; presence of large obstructing distal ureteral stone.

Potential Complications

Potential complications Include but not limited to Edema, irritation, inflammation & ureteral avulsion.

Warnings And Precautions

- Manipulation of device requires Ultrasound / Fluoroscopy / Other image guiding techniques.
- Do not use if product is damaged or expired and reuse of this device leads to potential complications.
- This device is intended for use by physicians trained and experienced in diagnostic and endoscopic techniques, standard techniques for placement of endoscopic devices should be employed.
- On completion of the procedure, dispose off the device and its components, if any, as per institutional and local regulatory guidelines for bio-hazardous medical waste disposal.
-
- Not intended for permanent use.
- Never use undue force for placement of this device.
- Use only compatible endoscopic devices.
- Do not re-sterilize the device. It might show loss in structural integrity.
- Refer to the device label, for technical specifications.
- Obturator shall be locked while inserting the sheath through the ureter.
- Do not force the device when resistance occurs while usage.
- Device is intended for single use.

Instructions For Use

Conventional Placement Technique

- Prior to placement, activate the hydrophilic coating by immersing all components in sterile water or isotonic solution.
- This will make the sheath lubricious.
- Place 0.038inch guidewire of 150 cm length into the ureter to establish a working tract using standard procedure.
- Ensure the obturator is securely locked on to the Sheath adapter, allowing the obturator and sheath assembly to be placed as a single unit.
- Advance the distal tip of the obturator and sheath assembly over the guidewire into the ureter.
- Confirm proper placement using imaging guidance.
- While maintaining sheath position, unlock the obturator hub and remove the obturator. DO NOT advance the sheath without the obturator in place.
- The desired endoscope or instruments can now be used through the ureteral access sheath as needed.
- On completion of the access procedure, gently withdraw the device.

RAPID RELEASE TECHNIQUE

- Prior to placement, activate the hydrophilic coating by removing the dilator from the access sheath and immersing all components in sterile water or isotonic solution. This will make the sheath lubricious.
- Place an 0.038" inch guidewire of 150 cm length into the ureter to establish a working tract using standard procedure.
- Ensure, the dilator is securely locked on to the Sheath adapter, allowing the dilator or sheath assembly to be placed as a single unit.
- Advance the dilator and access sheath assembly over the guide wire from distal end of the dilator to the floppy cut of the dilator (Fig.1)
- Grasp the sheath just below the instrument adapter and advance the distal tip of the dilator/sheath assembly over the wire guide and into the ureter.
- Confirm the dilator/sheath assembly is properly placed via fluoroscopy.
- While holding the access sheath in position, unlock the fitting and remove the dilator in one slow, continuous pull. (Fig. 2) Removal of the dilator will release the wire guide via the slit on the distal tip. The wire guide will remain outside the working channel to maintain access. (Fig. 3)
- Introduce the desired instrument as needed.
- On completion of the access procedure, gently withdraw the device. Confirm the dilator/sheath assembly is properly placed via fluoroscopy.
- While holding the access sheath in position, unlock the fitting and remove the dilator in one slow, continuous pull. (Fig. 2) Removal of the dilator will release the wire guide via the slit on the distal tip. The wire guide will remain outside the working channel to maintain access. (Fig. 3)
- Introduce the desired instrument as needed.
- On completion of the access procedure, gently withdraw the device.

Fig-1


Fig-2

Fig-3

Target Population Group

Adult Male and Female population groups.

Supply and Storage

The device is packed in primary pack made of polyethylene and the secondary package made of polyethylene and Tyvek. Supplied sterilized by Ethylene Oxide gas in peel open package. Device shall be stored at temperature range from 20°C to 40°C and away from direct sunlight.

Report Any Serious Incident

If problems occur using this device or any serious incidents or adverse events, please contact our customer relations department or call to the contact provided in our manufacturing address and competent authority of the Member State in which the user and/or patient is established.

Symbol Glossory And Definations

Standard Reference:

ISO 15223-1 Medical devices Symbols to be used with information to be supplied by the manufacturer

Symbol	Symbol Title
	Do not re-sterilize
	Do not use if package is damaged and consult instructions for use
	Consult instructions for use
	Temperature limit
	Keep away from sunlight
	Do not re-use
	Not made with natural rubber latex
	Sterilized using ethylene oxide
	Non-pyrogenic
	Keep dry
	Medical device
	Single sterile barrier system with protective packaging inside
	Does not contain DEHP
	Unique device identifier
	Date of Manufacture
	Use-by date
	Batch code

	Catalogue number
	Authorized Representative in the European Community/ European Union
	Manufacturer
	Caution
	Prescription

-  only
- This device is to be sold by/on the directions of a physician. Sterile (ETO) if package is unopened or undamaged.
 - Read all the instructions before use.
 - On completion of the procedure, dispose off device and its components if any as per institutional and local regulatory guidelines for bio-hazardous medical waste dispose.

Authorized Representative in the European Community
Arazy Group (Ireland) Ltd
19 Baggot Street Lower, Dublin 2,
D02 X658, Ireland
E-mail: Ireland@arazygroup.com

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Blue Neem Medical Devices Pvt Ltd