

Grasping forceps

Size (Fr)	Length (cm)	Model	Type
2.5, 3, 4, 5, 6, 7, 8, 9 & 12	37, 40, 45, 60, 70, 110, 115 & 120	Tri prong Biprong Alligator	Rigid-URS, Rigid-PCNL, Rigid-Cystoscope, Flexible

Intended Purpose & Duration

Used to remove the stones and retrieval of foreign object under direct vision up to 60 minutes.

Contraindications

There are no known absolute contraindications for this device. However, the following warning is listed as a relative contraindication;

- Some stones may be too large to remove endoscopically with a forceps.
- Do not use the forceps, if stone is too large to be removed endoscopically or to be held by the forceps.
- To avoid stone impaction, fluoroscopy or x-ray should be used to determine the size of the stone.

Potential Complications

Complications may occur at any time during or following the procedure. Possible complications include: Edema, Urosepsis, Evulsion, Perforation of the ureter & Iatrogenic related infections.

Warnings And Precautions

- The device is conductive. Avoid contact with electrified instruments and generators.
- Do not reuse the device. Reuse can cause severe infection to the patient. Intended use may fail on reuse. Do not use if package is opened or found defected. Do not use the product if there is suspect that the device is not sterile.
- Device is intended to be used by physicians trained and experienced in the endoscopic urological procedures with support from qualified technicians using aseptic technique. On completion of the procedure, dispose off the device and its components, if any, as per institutional and local regulatory guidelines for bio-hazardous medical waste disposal.
- Device is not intended for permanent use.

Instructions For Use

1. Choose a forcep of appropriate size, length and type depending on the procedure and size, location of the stone or foreign body [Object].
2. Through the working channel of the endoscope, pass the forcep to the site of the Object.
3. When the Object position is reached, press or operate the handle to open the forcep, catch the Object and release or operate the handle to secure the Object.
4. Carefully retrieve the object through the scope or along with the scope.

Intended population group

Patients who require removal of foreign object during endoscopic procedure.

Supply and Storage

The device is packed in package made of polyethylene and Tyvek for device up to length 60 cm. More than 60 cm length device is packaged in tray and packaged in Tyvek pouch. Supplied sterilized by Ethylene Oxide gas in peel open package. Device shall be stored at temperature range from 20°C to 40°C and away from direct sunlight.

Report Any Serious Incident

If problems occur using this device or any serious incidents or adverse events, please contact our customer relations department or call to the contact provided in our manufacturing address and competent authority of the Member State in which the user and/or patient is established.

Symbol Glossory And Definations

Standard Reference:

ISO 15223-1 Medical devices — Symbols to be used with information to be supplied by the manufacturer

Symbol	Symbol Title
	Do not re sterilize
	Do not use if package is damaged and consult instructions for use
	Consult instructions for use
	Temperature limit
	Keep away from sunlight
	Do not re-use
	Not made with natural rubber latex
	Sterilized using ethylene oxide
	Non-pyrogenic
	Keep dry
	Medical device
	Single sterile barrier system with protective packaging inside
	Does not contain DEHP
	Unique device identifier
	Date of Manufacture
	Use-by date
	Batch code

	Catalogue number
	Authorized Representative in the European Community/ European Union
	Manufacturer
	Caution
	Prescription



- This device is to be sold by/on the directions of a physician. Sterile (ETO) if package is unopened or undamaged.
- Read all the instructions before use.
- On completion of the procedure, dispose off device and its components if any as per institutional and local regulatory guidelines for bio-hazardous medical waste dispose.



Authorized Representative in the European Community
 Arazy Group (Ireland) Ltd
 19 Baggot Street Lower,
 Dublin 2, D02 X658, Ireland
 E-mail: Ireland@arazygroup.com



Manufacturer
Blue Neem Medical Devices Pvt Ltd
 Plot Nos 270 & 271, Road No 5, Harohalli Industrial Area
 II Phase, Kanakapura Taluk, Ramanagara
 Karnataka - 562112, India.
 Phone: +91 80 29761336
 E-Mail: contact@blueneem.com
 Website: www.blueneem.com