Hydro Twister[®] Guidewire

Device Description

Hydrotwister guidewires are made of a Nitinol core mandril inside and outer Polymer sleeve. Hydrotwister guidewire is thin wire with flexible hydrophilic tip on proximal end and rigid hydrophilic tip available in different diameters, lengths and straight flexible tip configurations that can be inserted into a tract to facilitate the placement of a medical device or instrument such as an endoscope, catheter, stent or dilator.

The surface of the wire is made in such a way to facilitate coaxial use with other devices. Please refer to the product label for product specifications.

Hydrotwister Guidewire is available in following tip configuration:

- Straight tip

Size (inch)	Length (cm)	Variant
0.025, 0.032, 0.035 & 0.038	150 & 260	Hydrotwister Guidewire Hydrotwister Plus Guidewire

Intended Purpose

Guidewire is intended to gain access to desired anatomical location and to assist in placement of medical devices such as catheters, stents and dilators during diagnostic or interventional procedures.

Indications For Use & Contraindications

Indicated to be used to gain the access to desired anatomical site to assist in the placement of medical device. No known contraindications.

Potential Complications

Potential complications include but not limited to Tissue perforation, Infection, Tissue trauma & Edema.

Warnings And Precautions

- Single use only. Do not reuse, reuse can cause severe infection. Intended use may fail on reuse and might cause kink of wire on reuse.
- Use on or before the shelf life printed on label.
- Do not resterilize & use only if packaging is intact.
- This device is not intended to be used for coronary arteries and other delicate tortuous anatomical regions.
- Avoid forceful angulation. Avoid manipulating or withdrawing the guidewire back through metal cannula or needle. A sharp edge may scrape or shear material from wire.
- Do not alter the tip configuration manually. It may damage the guidewire.
- This device is intended to be used by trained physicians or trained technicians.
- Use medical imaging when manipulating the guidewire. Do not advance or manipulate the guidewire without the visual evidence of the corresponding movement of the distal tip.
- When guidewire is used in association with another device, compatibility shall be checked before using in the procedure.
- Use only flexible end of the guidewire as initial inserted end. Inserting from rigid end may cause the damage to tissue or perforation.
- Inspect the wire for kinks and other damages on sleeve prior to use.
- Device is intended for transient and single use.

Instructions For Use

- Open the device package and remove the guidewire using aseptic technique.
- Locate the flexible tip of the guidewire and ensure the guidewire compatibility with other device, then insert or advance the wire through the device or orifice and position the flexible tip of the guidewire in desired anatomical location.
- Use only medical imaging techniques for manipulating the guidewire.
- Perform the procedure and retrieve the guidewire gently.
- Discard the wire after the procedure.

Supply and Storage

The device is assembled in protective dispenser tube and the secondary package made of polyethylene and Tyvek. Supplied sterilized by Ethylene Oxide gas in peel open package. Device shall be stored at temperature range from 20°C to 40°C and away from direct sunlight.

Report Any Serious Incident

If problems occur using this device or any serious incidents or adverse events, please contact our customer relations department or call to the contact provided in our manufacturing address and competent authority of the Member State in which the user and/or patient is established.

Symbol Glossory And Definations

Standard Reference:

ISO 15223-1 Medical devices — Symbols to be used with information to be supplied by the manufacturer

Symbol	Symbol Title
	Do not re-sterilize
	Do not use if package is damaged and consult instructions for use
	Consult instructions for use
	Temperature limit
	Keep away from sunlight
	Do not re-use
	Not made with natural rubber latex
	Sterilized using ethylene oxide
	Non-pyrogenic
	Keep dry
	Medical device
	Single sterile barriersystem with protective packaging inside
	Does not contain DEHP
	Unique device identifier
	Date of Manufacture
	Use-by date
	Batch code

	Catalogue number
	Authorized Representative in the European Community/ European Union
	Manufacturer
	Caution
	Prescription

- This device is to be sold by/on the directions of a physician. Sterile (ETO) if package is unopened or undamaged.
- Read all the instructions before use.
- On completion of the procedure, dispose off device and its components if any as per institutional and local regulatory guidelines for bio-hazardous medical waste dispose.

Authorized Representative in the European Community
Arazy Group (Ireland) Ltd
19 Baggot Street Lower,
Dublin 2, D02 X658, Ireland
E-mail: Ireland@arazygroup.com

Manufacturer
Blue Neem Medical Devices Pvt Ltd
Plot Nos 270 & 271, Road No 5, Harohalli Industrial Area
II Phase, Kanakapura Taluk, Ramanagara
Karnataka - 562112, India.
Phone: +91 80 29761336
E-Mail: contact@blueneem.com
Website: www.blueneem.com

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