Mono J Catheter

Device Description

Mono J Catheter is a flexible tube-like device made from polyurethane for introduction into urinary tract through an endoscope. Device is available in variety of French sizes, pigtail distal tip configurations for retention. Choice of device is based on clinical situation. Refer to label for technical specifications.

Device package includes (contents may varies):

- Radiopaque pigtail catheter
- Urine bag connector
- PTFE Wire guide (available only in catheter sets)

Size (Fr)	Length (cm)	Variant
4.7, 5, 6, 7 & 8	70	Open end & Closed end

Intended Purpose

Used for temporary drainage and irrigation of urine from the kidney for a period of up to 48 hours.

Indications For Use & Contraindications

Contraindicated for use in patients unsuitable for catheterization.

Potential Complications

Potential complications Include but not limited to haemorrhage, Edema, perforation of the kidney, ureter avulsion renal pelvis, ureter and/or bladder, urinary tract infection.

Warnings And Precautions

- Manipulation of device requires Ultrasound / Fluoroscopy / other image guiding techniques.
- Do not use if product is damaged or expired, reuse and resterilization of this device leads to potential complications and loss in structural integrity.
- This device is intended for use by physicians trained and experienced in diagnostic and endoscopic techniques, standard techniques for placement of endoscopic devices should be employed.
- On completion of the procedure, dispose off the device and its components, if any, as per institutional and local regulatory guidelines for bio-hazardous medical waste disposal.
- Not intended for permanent use.
- Use only compatible endoscopic devices.
- Avoid bending or kinking of catheter prior to placement. Doing so could damage the integrity of the catheter and result in patient injury.

Instructions For Use

Procedure For Closed End

- Insert the compatible guidewire with stiff tip into the straight end of the catheter and stretch the catheter pigtail coil at the other end.
- Advance the catheter and guidewire until the tip reaches the renal pelvis.
- Once correct placement within the renal pelvis is confirmed by fluoroscopic guidance, remove the guidewire carefully.

Procedure For Open End

- Insert soft tip of the compatible guidewire into endoscope, access ureteral orifice and pass the guidewire up to the renal Pelvis.
- Pass the catheter over guidewire through the endoscope maintaining guidewire position.
- Push the catheter until it reaches renal pelvis under fluoroscopic guidance.
- Once correct placement within the renal pelvis is confirmed by fluoroscopic guidance, remove the guidewire carefully.
- Remove the endoscope maintaining the position of the catheter.

NOTE: This Device is detectable by X-Ray or other means of fluoroscopy.

Catheter Removal:

Before the intended duration completes, gently withdraw the catheter manually either by hand or help of forceps.

Product Recommendations:

Recommended wire guide for the usage of this product is as follows:

4.7Fr: 0.032inch; 05Fr: 0.035inch; 06, 07 & 08Fr: 0.038inch.

Supply and Storage

The device is packed in primary pack made of polyethylene and the secondary package made of polyethylene and Tyvek. Supplied sterilized by Ethylene Oxide gas in peel open package. Device shall be stored at temperature range from 20°C to 40°C and away from direct sunlight.

Report Any Serious Incident

If problems occur using this device or any serious incidents or adverse events, please contact our customer relations department or call to the contact provided in our manufacturing address and competent authority of the Member State in which the user and/or patient is established.

Symbol Glossory And Definations

Standard Reference:

ISO 15223-1 Medical devices — Symbols to be used with information to be supplied by the manufacturer

Symbol	Symbol Title
	Do not re-sterilize
	Do not use if package is damaged and consult instructions for use
	Consult instructions for use
	Temperature limit
	Keep away from sunlight
	Do not re-use
	Not made with natural rubber latex
	Sterilized using ethylene oxide
	Non-pyrogenic
	Keep dry
	Medical device
	Single sterile barrier system with protective packaging inside
	Does not contain DEHP
	Unique device identifier
	Date of Manufacture
	Use-by date
	Batch code

	Catalogue number
	Authorized Representative in the European Community/ European Union
	Manufacturer
	Caution
	Prescription

- This device is to be sold by/on the directions of a physician. Sterile (ETO) if package is unopened or undamaged.
- Read all the instructions before use.
- On completion of the procedure, dispose off device and its components if any as per institutional and local regulatory guidelines for bio-hazardous medical waste dispose.

Authorized Representative in the European Community
Arazy Group (Ireland) Ltd
19 Baggot Street Lower,
Dublin 2, D02 X658, Ireland
E-mail: Ireland@arazygroup.com

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Blue Neem Medical Devices Pvt Ltd