N Quadra[®], N Quadra[®] Z

N Quadra[®] Helical

N Quadra[®] Helical Z

N Quadra[®] Flexible

Device Description

N Quadra[®], N Quadra[®] Z, N Quadra[®] Helical, N Quadra[®] Helical Z, N Quadra[®] Flexible are four wire nitinol stone baskets available with adjustable handle which allows controlled assured basket expansion and contraction in a range of French sizes, lengths and basket diameters.

Blueneem Baskets are available in flexible and semi rigid working length configurations with variable basket diameters . Please refer to the product label for product specifications.

Model	Size (Fr)	Length (cm)
N Quadra [®]	2.5, 3, 4 & 5	70,90 & 120
N Quadra [®] Z	2.5, 3, 4 & 5	70,90 & 120
N Quadra [®] Helical	2.5, 3, 4 & 5	70,90 & 120
N Quadra [®] Helical Z	2.5, 3, 4 & 5	70,90 & 120
N Quadra [®] Flexible	1.8, 1.9, 2, 2.2 & 2.4	115 & 120

Intended Purpose

Stone Retrieval Basket is intended to be used in an endoscopic procedure for removal of urinary stones or foreign objects from the Ureter or the Kidney.

Indications For Use & Contraindications

Indicated to be used to retrieve urinary stones or foreign objects form ureter or kidney. There are no known absolute contraindications but following is listed but not limited.

Do not use the basket to grip the stone when its size is too large to remove endoscopically.

To avoid stone impaction fluoroscopy or x-ray should be used to determine the size of the stone.

Potential Complications

Complications may occur at any time during or following the procedure. Possible complications include: Oedema, Urosepsis, Perforation of the ureter & Iatrogenic related infections.

Warnings And Precautions


- Manipulation of the device needs ultrasound or fluoroscopic guidance.
- If stone size is too large, break the stone into fragments and then retrieve, otherwise it may cause damage to the basket.
- To avoid stone impaction, fluoroscopy or X-Ray should be used to assess the size of the stone.
- The device is conductive. Avoid contact with electrified instrument.
- Do not use if product is damaged or expired or package integrity is compromised.
- Do not reuse the device. Re use/ Resterlization of the device can compromise structural integrity of the device thereby causing serious harm to the patient. Reuse can contaminate the device causing serious harm to the patient by way of infection.
- Do not open the basket or operate the device handle when the tip of the opening side of the basket is still inside the Endoscope working channel. It may lead to failure of the device.
- For semi rigid baskets do not bend the working length of the basket. As it may lead to kink or break the device.

Precautions

- Device use is intended to be carried out by trained physicians and technicians who were trained in endoscopic urological procedures. Standard endoscopic procedures should be employed.
- Do not use excessive force to manipulate the device, it may cause damage to the device.
- Only intended for single use.

Instruction of Use

- Remove the stone basket from outer package and move the handle position to “Close” position to close the basket as shown in below figures.
- Insert the device into the scope and under direct visualization, locate the stone and advance the basket to the location just proximal to the stone.
- Slide the handle lever to open position as shown in below figures for different types of handles to deploy the basket near the stone.
- Maneuver the device until the stone is within the basket
- Once correctly positioned, draw the basket down around the stone and operate the handle as shown in below figures.
- Maintaining position of the stone within the basket, retract the scope and the basket

Target Population Group

Adult male and female individuals.

Supply and Storage

The device is packed in blister tray sealed with medical grade Tyvek pouch. Supplied sterilized by Ethylene Oxide gas in peel open package. Device shall be stored at temperature range from 20°C to 40°C and away from direct sunlight.

Report Any Serious Incident

If problems occur using this device or any serious incidents or adverse events, please contact our customer relations department or call to the contact provided in our manufacturing address and competent authority of the Member State in which the user and/or patient is established.

Symbol Glossary And Definitions

Standard Reference:

ISO 15223-1 Medical devices Symbols to be used with information to be supplied by the manufacturer

Symbol	Symbol Title
	Do not re-sterilize
	Do not use if package is damaged and consult instructions for use
	Consult instructions for use
	Temperature limit
	Keep away from sunlight
	Do not re-use
	Not made with natural rubber latex
	Sterilized using ethylene oxide
	Non-pyrogenic
	Keep dry
	Medical device
	Single sterile barrier system with protective packaging inside
	Does not contain DEHP
	Unique device identifier
	Date of Manufacture
	Use-by date
	Batch code

	Catalogue number
	Authorized Representative in the European Community/ European Union
	Manufacturer
	Caution
	Prescription

- This device is to be sold by/on the directions of a physician. Sterile (ETO) if package is unopened or undamaged.
- Read all the instructions before use.
- On completion of the procedure, dispose off device and its components if any as per institutional and local regulatory guidelines for bio-hazardous medical waste dispose.

Authorized Representative in the European Community
Arazy Group (Ireland) Ltd
19 Baggot Street Lower,
Dublin 2, D02 X658, Ireland
E-mail: Ireland@arazygroup.com

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Blue Neem Medical Devices Pvt Ltd