Nephrostomy Balloon Dilator

Device Description

Device package contains (contents may vary)

- Balloon dilator with syringe
- Balloon dilator with inflation device
- Balloon dilator with syringe and Amplatz sheath
- Balloon dilator with inflation device and Amplatz sheath

Device contents may vary based on requirements.

Intended Purpose

Intended for fascial dilation, intended for stretching or enlarging the fascial tissue percutaneously prior to performing nephrostomy procedure.

Indications For Use & Contraindications

This device is contraindicated in the presence of conditions which create unacceptable risk during nephrostomy procedure.

Potential Complications

Potential complications include but not limited to tissue trauma, bleeding, perforation, infection, urine leakage & pain.

Warnings And Precautions

- This device is intended for use by physicians trained and experienced in techniques for balloon dilation or catheterization and this device should not be used without comprehensive knowledge of the indications, techniques and risks of the procedure.
- Always inflate the balloon with sterile liquid. Do not inflate the balloon with any kind of air or carbon dioxide or other gas. Do not over inflate or exceed maximum rated burst pressure, as excessive pressure may cause the balloon to rupture.
- Refer to product label for appropriate inflation details.
- If loss of pressure within the balloon occurs during inflation or if balloon ruptures during dilatation, immediately discontinue the procedure. Deflate the balloon carefully. Do not attempt to reinflate the balloon. Carefully withdraw the device if failure found during procedure.
- Do not exceed the maximum rated burst pressure (listed on label) for this balloon device. Refer label for device specifications.
- Use only compatible wire guide as specified on the device label.
- Do not use if product is damaged or expired, reuse and re-sterilization of this device leads to loss in structural integrity and performance. Device is intended for single use.
- On completion of the procedure, dispose of the device as per local regulatory guidelines for bio-hazardous medical waste disposal.

Instruction of Use

Procedure For Balloon Dilator Without Amplatz Sheath

- Perform initial access with puncture needle for the purpose of sequential dilation using standard procedure.
- Remove the sleeve from the device and discard.
- Lock the syringe on to the side arm stop cock of the balloon catheter and aspirate until all air is removed and close the stop cock.
- Using imaging guidance, pass a guidewire of 0.035 inch to the desired distance of the tract to be made.
- Pass the balloon dilator over the guidewire and position the balloon with the help of radiopaque markers under imaging guidance, ensure proper placement of the balloon at the site of dilation.
- Fill the syringe or use the inflation device if provided, with contrast medium and eliminate air from the syringe or inflation device.
- Open the stop cock on the dilator, connect the syringe or inflation device and inflate the balloon with the contrast medium. CAUTION: Inflate the balloon only with sterile medium. Do not inflate with air or any gas. To ensure proper regulation of balloon pressure, use of balloon inflation device with pressure gauge is recommended. Do not exceed the recommended maximum pressure mentioned on the product label.

- On reaching desired balloon pressure, close the stop cock and leave the inflated balloon at place for at least 30 seconds to allow adequate dilation.
- Once desired dilation is achieved, open the stop cock and aspirate the liquid from the balloon using the syringe or inflation device.
- Gently remove the dilator from the nephrostomy tract.

Procedure For Balloon Dilator With Amplatz Sheath

- Remove the sleeve from the device and discard.
- Lock the syringe on to the side arm stop cock of the balloon catheter and aspirate until all air is removed and close the
- Using fluoroscopic guidance, pass a guidewire of 0.035Inch to the desired distance of the tract to made.
- Introduce the balloon dilator through the lumen of the amplatz sheath.
- Pass the balloon dilator over the guidewire and position the balloon with the help of radiopaque markers under fluoroscopic guidance, Please Ensure proper placement of the balloon at the site of dilation.
- Fill the syringe or use the inflation device if provided, with contrast medium and eliminate air from the syringe or inflation device.
- Open the stop cock on the dilator, connect the syringe or inflation device and inflate the balloon with the contrast medium. CAUTION: Inflate the balloon only with sterile medium. Do not inflate with air or any gas. To ensure proper regulation of balloon pressure, use of balloon inflation device with pressure gage is recommended. Do not exceed the recommended max pressure mentioned on the product label.
- On reaching desired balloon pressure, close the stop cock and leave the inflated balloon at place for at least 30 seconds to allow adequate dilation.
- Once desired dilation is achieved, advance the previously loaded amplatz sheath onto the balloon to create the nephrostomy tract, after that open the stop cock and aspirate the liquid from the balloon using the syringe or inflation device.
- Gently remove the dilator from the nephrostomy tract.

Target Population Group

Adult Male and Female population groups.

Supply and Storage

The device is assembled in protective dispenser tube and the secondary package made of polyethylene and Tyvek. Supplied sterilized by Ethylene Oxide gas in peel open package. Device shall be stored at temperature range from 20°C to 40°C and away from direct sunlight.

Report Any Serious Incident

If problems occur using this device or any serious incidents or adverse events, please contact our customer relations department or call to the contact provided in our manufacturing address and competent authority of the Member State in which the user and/or patient is established.

Symbol Glossory And Definations

Standard Reference:

ISO 15223-1 Medical devices - Symbols to be used with information to be supplied by the manufacturer

Symbol	Symbol Title
	Do not re-sterilize
	Do not use if package is damaged and consult instructions for use
	Consult instructions for use
	Temperature limit
	Keep away from sunlight
	Do not re-use
	Not made with natural rubber latex
	Sterilized using ethylene oxide
	Non-pyrogenic
	Keep dry
	Medical device
	Single sterile barrier system with protective packaging inside
	Does not contain DEHP
	Unique device identifier
	Date of Manufacture
	Use-by date
	Batch code

	Catalogue number
	Authorized Representative in the European Community/ European Union
	Manufacturer
	Caution
	Prescription

- This device is to be sold by/on the directions of a physician. Sterile (ETO) if package is unopened or undamaged.
- Read all the instructions before use.
- On completion of the procedure, dispose off device and its components if any as per institutional and local regulatory guidelines for bio-hazardous medical waste dispose.

Authorized Representative in the European Community
Arazy Group (Ireland) Ltd
19 Baggot Street Lower, Dublin 2,
D02 X658, Ireland
E-mail: Ireland@arazygroup.com

Manufacturer
Blue Neem Medical Devices Pvt Ltd
Plot Nos 270 & 271, Road No 5, Harohalli Industrial Area
II Phase, Kanakapura Taluk, Ramanagara
Karnataka - 562112, India.
Phone: +91 80 29761336
E-Mail: contact@blueneem.com
Website: www.blueneem.com

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