PCN Catheter

Device Description

PCN Catheter is flexible tube made of polyurethane, intended to be inserted into but not limited into the renal pelvis percutaneously and to gain access into the upper urinary tract. The tip of the device has a pigtail shape for retention. The retention means have eye type ports to drain the fluids. Device packaging is provided with plastic/metal retention loop cannulae & urine bag connector. For product specifications refer to the label. PCN Catheters package includes (set components may vary):

PCN Catheters (Standard)

package includes (set components may vary)

- Radiopaque pigtail catheter (with or without Hydrophilic coating)
- Plastic/metal cannulae
- Plastic sheath
- Urine bag connector

PCN Catheters with locking

mechanism package includes (set components may vary)

- Radiopaque pigtail catheter with locking mechanism (with or without Hydrophilic coating)
- Plastic/ metal cannulae
- Plastic sheath
- Urine bag connector

Size (Fr)	Length (cm)	Model
6, 7, 8, 8.5, 9, 10, 12, 14, 16 & 18	20, 22, 25, 30, 35, 40 & 45	PCN Catheter (with or without Locking Mechanism). Available with or without Hydrophilic coating.

Intended Purpose & Duration

PCN catheters are intended to be inserted percutaneously into the renal pelvis and to gain access to the upper urinary tract to perform external drainage of urine from the renal pelvis by seldinger access technique for a short-term period of up to 30 days.

Indications For Use & Contraindications

When diagnosed to place drainage catheter percutaneously in kidney for renal drainage. No known contraindications. Device should not be used in the conditions of Bleeding Diathesis and uncontrolled hypertension, Anticoagulant use.

Potential Complications

Bleeding (including chronic hematuria, post-procedural, intraprocedural, puncture site), Catheter dislocation/dislodgement, Catheter malposition/transcolonic catheter placement, Catheter occlusion/obstruction, Damage to adjacent organs, Urinary tract infection, Mechanical catheter damage and pain.

Warnings And Precautions

- If a catheter has become malpositioned or if drainage ceases, the catheter should be promptly exchanged or removed.
- The split away Plastic Sheath, if retained, is not be used as introducer sheath for any purpose.
- Do not use if product is damaged or expired and do not reuse the device.
- Use only accessories supplied along with package.
- Manipulation of the product needs ultrasound or fluoroscopic guidance.
- If the device is hydrophilic coated, activate the coat by using sterile water or saline or other isotonic solution to make the catheter insertion easier.
- Dispose off the device and its components, if any, as per institutional and local regulatory guidelines for bio-hazardous medical waste disposal.
- This device is intended for use by physicians trained and experienced in diagnostic and interventional techniques, standard techniques for placement of percutaneous drainage catheters should be employed.

- Periodic evaluation of the catheter is advised; catheter should not remain indwelled for more than 30 days after deployment.
- Do not forcefully advance any components during removal or placement of the catheter. Carefully remove the components if any resistance is occurred.
- Reuse and resterilization of the device is strictly prohibited.

Catheter Placement Procedure


For **Catheters supplied with hydrophilic coating**: Activate the hydrophilic coating, by wetting the catheter with sterile water or saline. For best results, keep catheter surface wet during placement.

Catheter (without locking mechanism)

1. Perform procedure using standard procedures of puncture, compatible guidewire placement and dilatation.
2. Use the peel away plastic sheath provided on the catheter to straighten the pigtail of the catheter and then insert the plastic or metal cannulated stiffener provided and twist the hub to lock. Then remove the peel away plastic sheath.
3. Load the pigtail end of the catheter over the guidewire.
4. Gradually advance the pigtail end well into the collecting system and remove the stiffener.
5. Confirm the pigtail position using fluoroscopic guidance.
6. Holding the shaft of the catheter securely in position remove the guidewire carefully.
7. Secure the catheter to the catheter-fixing disc or to the skin.
8. Using the Urine bag connector, connect the catheter to the urine collection bag.

Catheter (with locking mechanism)

1. Perform procedure using standard procedures of puncture, compatible guidewire placement and dilatation.
2. Use the peel away plastic sheath provided with the catheter to straighten the pigtail of the catheter.
3. Then insert the stiffener provided into the catheter and twist to lock with catheter hub.
4. Remove the Peel away plastic Sheath on the catheter. When inserting the stiffener into the catheter with locking mechanism, hold the string during stiffener insertion to avoid bunching or tangling of suture. Luer lock should not be tighten along with stiffener, to do so may cause kinking of catheter.
5. Then pass the pigtail end of the catheter over the external end of the guidewire.
6. Gradually advance the pigtail end well into the collecting system and remove the stiffener.
7. Confirm the pigtail position using fluoroscopic guidance.
8. Holding the shaft of the catheter securely in position remove the guidewire carefully, confirm placement via free flow of urine. On confirmation, remove the stiffener. The catheter will now acquire pigtail shape.
9. After that pull the locking string as shown in image A.
10. After pulling the string, wind the string firmly around the locking hub as shown in image B.
11. Then lock the red clip to the locking hub, therefore the pigtail locks firmly inside the body as shown in image C.
12. Traction on the locking string, if present, should be sufficient to ensure adequate retention of the tip, but should not be overly tight. Secure the catheter to the catheter-fixing disc or to the skin.
13. Using the Urine bag connector, connect the catheter to the urine collection bag.


Catheter Removal For Standard Type

- Detach the catheter from catheter fixing disc or skin.
- Using the plastic stiffener supplied along with the catheter make the pigtail of catheter straight and remove from body carefully.
- If necessary, place a guidewire, into the collecting system to maintain access while removing the catheter.
- Remove the catheter gently and carefully.

Catheter Removal For Locking Mechanism Type

- For removal of Catheter from the body unlock the red clip from the locking hub of the catheter.
- Unwind the string in anti-clockwise, after that turn the locking hub in anti-clockwise direction to release the locking of the string.
- Using the plastic stiffener supplied along with the catheter make the pigtail of catheter straight and remove from body carefully.
- If necessary, place a guidewire, into the collecting system to maintain access while removing the catheter.

Product Recommendations

Stent Size (Fr)	Accepted Guidewire (inch)
6,7,8 & 8.5	0.035
9, 10, 12, 14, 16 & 18	0.038

Supply and Storage

The device is packed in primary pack made of polyethylene and the secondary package made of polyethylene and Tyvek. Supplied sterilized by Ethylene Oxide gas in peel open package. Device shall be stored at temperature range from 20°C to 40°C and away from direct sunlight.

Target Population Group

In patients who have, or might have, a urinary obstruction due to various reasons, an external drainage tube called a 'catheter' is commonly placed in the kidney in order to prevent or to temporarily relieve the obstruction.

Report Any Serious Incident

If problems occur using this device or any serious incidents or adverse events, please contact our customer relations department or call to the contact provided in our manufacturing address and competent authority of the Member State in which the user and/or patient is established.

Symbol Glossory And Definations

Standard Reference:

ISO 15223-1 Medical devices — Symbols to be used with information to be supplied by the manufacturer

Symbol	Symbol Title
	Do not re-sterilize
	Do not use if package is damaged and consult instructions for use
	Consult instructions for use
	Temperature limit
	Keep away from sunlight
	Do not re-use
	Not made with natural rubber latex
	Sterilized using ethylene oxide
	Non-pyrogenic
	Keep dry
	Medical device
	Single sterile barrier system with protective packaging inside
	Does not contain DEHP
	Unique device identifier
	Date of Manufacture
	Use-by date
	Batch code

	Catalogue number
	Authorized Representative in the European Community/ European Union
	Manufacturer
	Caution
	Prescription

- This device is to be sold by/on the directions of a physician. Sterile (ETO) if package is unopened or undamaged.
- Read all the instructions before use.
- On completion of the procedure, dispose off device and its components if any as per institutional and local regulatory guidelines for bio-hazardous medical waste dispose.

Authorized Representative in the European Community
Arazy Group (Ireland) Ltd
19 Baggot Street Lower,
Dublin 2, D02 X658, Ireland
E-mail: Ireland@arazygroup.com

Manufacturer
Blue Neem Medical Devices Pvt Ltd
Plot Nos 270 & 271, Road No 5, Harohalli Industrial Area
II Phase, Kanakapura Taluk, Ramanagara
Karnataka - 562112, India.
Phone: +91 80 29761336
E-Mail: contact@blueneem.com
Website: www.blueneem.com

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