PCNL Basket

Device Description

The PCNL basket is a stone extracting device with a rigid shaft and an expandable basket of wire made of Nitinol and shaft made of stainless steel. It is designed to extract urinary stones from the calices of kidney in PCNL procedure. It is introduced through the working channel of the rigid nephroscope.

Please refer to the product label for product specifications.

Intended Purpose

PCNL Basket is intended to be used with rigid nephroscope for extraction of stones from calices.

Contraindications

None Known

Warnings

- Manipulation of the device needs ultrasound or fluoroscopic guidance.
- If stone size is too large, break the stone into fragments and then retrieve, otherwise it may cause damage to the basket.
- The device is conductive. Avoid contact with electrified instrument.
- Do not use if product is damaged or expired and packaging integrity is compromised
- Do not reuse the device. Re use/ Resterlization of the device can compromise structural integrity of the device thereby causing serious harm to the patient. Reuse can contaminate the device causing serious harm to the patient by way of infection.
- Device use is intended to be carried out by trained physicians and technicians who were trained in endoscopic urological procedures. Standard endoscopic procedures should be employed.
- Do not use excessive force to manipulate the device, it may cause damage to the device.
- Only intended for single time and transient use.

Precautions

- Device use is intended to be carried out by trained physicians and technicians who were trained in endoscopic urological procedures. Standard endoscopic procedures should be employed.
- Do not use excessive force to manipulate the device, it may cause damage to the device.
- Not intended for permanent use. Only intended for transient use.

Instruction for Use

- During PCNL procedures, after access, dilatation and lithotripsy baskets will be used for retrieval of stones from the system.
- The Basket should be passed through the working channel of scope.
- By pressing the U shape handle the basket will open and catch the object on releasing the handle.
- Carefully retrieve the object through the scope or along with the scope.

Target Population Group

Adult Male and Female population groups.

Supply and Storage

The device is assembled in protective tube and the secondary package made of polyethylene and Tyvek. Supplied sterilized by Ethylene Oxide gas in peel open package. Device shall be stored at temperature range from 20°C to 40°C and away from direct sunlight.

Report Any Serious Incident

If problems occur using this device or any serious incidents or adverse events, please contact our customer relations department or call to the contact provided in our manufacturing address and competent authority of the Member State in which the user and/or patient is established.

Symbol Glossory And Definations

Standard Reference:

ISO 15223-1 Medical devices — Symbols to be used with information to be supplied by the manufacturer

Symbol	Symbol Title
	Do not re-sterilize
	Do not use if package is damaged and consult instructions for use
	Consult instructions for use
	Temperature limit
	Keep away from sunlight
	Do not re-use
	Not made with natural rubber latex
	Sterilized using ethylene oxide
	Non-pyrogenic
	Keep dry
	Medical device
	Single sterile barrier system with protective packaging inside
	Does not contain DEHP
	Unique device identifier
	Date of Manufacture
	Use-by date
	Batch code

	Catalogue number
	Authorized Representative in the European Community/ European Union
	Manufacturer
	Caution
	Prescription

- This device is to be sold by/on the directions of a physician. Sterile (ETO) if package is unopened or undamaged.
- Read all the instructions before use.
- On completion of the procedure, dispose off device and its components if any as per institutional and local regulatory guidelines for bio-hazardous medical waste dispose.

Authorized Representative in the European Community
 Arazy Group (Ireland) Ltd
 19 Baggot Street Lower,
 Dublin 2, D02 X658, Ireland
 E-mail: Ireland@arazygroup.com

Manufacturer
Blue Neem Medical Devices Pvt Ltd
 Plot Nos 270 & 271, Road No 5, Harohalli Industrial Area
 II Phase, Kanakapura Taluk, Ramanagara
 Karnataka - 562112, India.
 Phone: +91 80 29761336
 E-Mail: contact@blueneem.com
 Website: www.blueneem.com

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