Percutaneous Nephrostomy Drainage Catheter Set - Malecot

Device Description

Percutaneous Nephrostomy Drainage Catheter - malecot is flexible tube made of polyurethane, intended to be inserted into the renal pelvis percutaneously and to gain access into the upper urinary tract. The tip of the device has wings shaped configuration which enhances retention. For product specifications refer to the label.

Percutaneous Nephrostomy Drainage Set – Malecot package includes (set components may vary):

- Radiopaque pigtail catheter (with or without Hydrophilic coating)
- Plastic/metal cannulated stiffener
- Plastic peel away sheath
- Fascial dilators
- PTFE Guidewire
- Initial puncture (IP) needle
- Catheter fixing disc
- Scalpel blade
- Urine bag connector

Size (Fr)	Length (cm)	Model
8, 10, 12, 14, 16, 18 & 20	20, 22, 25, 30 & 35	All models available with and without hydrophilic coating

Intended Purpose & Duration

Nephrostomy Drainage Catheter is intended to be inserted percutaneously into the renal pelvis and to gain access to the upper urinary tract to perform external drainage of urine from the renal pelvis for a short-term period of up to 30 days.

Indications For Use & Contraindications

When diagnosed to place drainage catheter percutaneously in kidney for renal drainage. No known contraindications. Device should not be used in the conditions of Bleeding Diathesis and uncontrolled hypertension, Anticoagulant use.

Potential Complications

Bleeding (including chronic hematuria, post-procedural, intraprocedural, puncture site), Catheter dislocation/dislodgement, Catheter malposition/transcolonic catheter placement, Catheter occlusion/obstruction, Damage to adjacent organs, Urinary tract infection, Mechanical catheter damage and pain.

Warnings And Precautions

- If a catheter has become malpositioned or if drainage ceases, the catheter should be promptly exchanged or removed.
- The split away Plastic Sheath, if retained, is not be used as introducer sheath any purpose.
- Do not use if product is damaged or expired and do not reuse the device.
- Use only accessories supplied along with package.
- Manipulation of the product needs ultrasound or fluoroscopic guidance.
- If the device is hydrophilic coated, activate the coat by using sterile water or saline or other isotonic solution to make the catheter insertion easier.
- Retain the plastic/metal cannulated straightener to straighten the catheter loop to remove the catheter.

- Dispose off the device and its components, if any, as per institutional and local regulatory guidelines for bio-hazardous medical waste disposal
- This device is intended for use by physicians trained and experienced in diagnostic and interventional techniques, standard techniques for placement of percutaneous drainage catheters should be employed.
- Periodic evaluation of the catheter is advised; catheter should not remain indwelled for more than 30 days after deployment.
- Do not forcefully advance any components during removal or placement of the catheter. Carefully remove the components if any resistance is occurred.
- Reuse and resterilization of the device is strictly prohibited.

Instructions For Use

For **Catheters supplied with Hydrophilic coating**: Activate the hydrophilic coating, by wetting the catheter with sterile water or saline. For best results, keep catheter surface wet during placement.

- By retrograde pyelogram, ultrasound or Computed Tomography(CT) scan, localize the kidney/collection system to be drained.
- After pyelogram, Identify and mark the skin site overlaying the kidney/collection system. Anesthetize the marked skin site.
- With the help of scalpel, make a small incision.
- Under fluoroscopic guidance, introduce the IP needle provided along with set, at right angles to the skin at the marked site. After achieving placement, remove the inner trocar of the needle. Confirm placement via free flow of urine.
- On confirmation, pass the guidewire through the outer sheath of the IP Needle.
- On reaching the desired position, remove the outer sheath of the needle while maintaining guidewire position.
- To facilitate passage of the catheter, dilate the musculofascial track using dilators in sequence from smallest to the largest.
- Use the stiffener provided to close the Malecot wings and twist the hub to lock. Then pass the wings end of the catheter over the external end of the guidewire. Gradually advance the wings end well into the collecting system. Confirm the position using fluoroscopic guidance.
- Unlock the stiffener with a gentle twist and remove the Straightner from the catheter to open the Malecot wings.
- Holding the shaft of the catheter securely in position remove the guidewire carefully.
- Tape or suture the catheter fixing disc to the skin.
- Secure the catheter to the catheter fixing disc.
- Using the urine bag connector, connect the catheter to the urine bag.

Catheter Removal Procedure

- Detach the catheter from catheter fixing disc.
- Gently pull the catheter carefully to remove.
- If necessary, place a guidewire, into the kidney/collecting system to maintain access.

Supply and Storage

The device is packed in such a way that, all the contents are assembled in tray and the secondary package made of polyethylene and Tyvek. Supplied sterilized by Ethylene Oxide gas in peel open package. Device shall be stored at temperature range from 20OC to 40OC and away from direct sunlight.

Target Population Group

In patients who have, or might have, a urinary obstruction due to various reasons, an external drainage tube called a 'catheter' is commonly placed in the kidney in order to prevent or to temporarily relieve the obstruction

Report Any Serious Incident

If problems occur using this device or any serious incidents or adverse events, please contact our customer relations department or call to the contact provided in our manufacturing address and competent authority of the Member State in which the user and/or patient is established.

Symbol Glossory And Definations

Standard Reference:

ISO 15223-1 Medical devices — Symbols to be used with information to be supplied by the manufacturer

Symbol	Symbol Title
	Do not re-sterilize
	Do not use if package is damaged and consult instructions for use
	Consult instructions for use
	Temperature limit
	Keep away from sunlight
	Do not re-use
	Not made with natural rubber latex
	Sterilized using ethylene oxide
	Non-pyrogenic
	Keep dry
	Medical device
	Single sterile barrier system with protective packaging inside
	Does not contain DEHP
	Unique device identifier
	Date of Manufacture
	Use-by date
	Batch code

	Catalogue number
	Authorized Representative in the European Community/ European Union
	Manufacturer
	Caution
	Prescription

- This device is to be sold by/on the directions of a physician. Sterile (ETO) if package is unopened or undamaged.
- Read all the instructions before use.
- On completion of the procedure, dispose off device and its components if any as per institutional and local regulatory guidelines for bio-hazardous medical waste dispose.

Authorized Representative in the European Community
Arazy Group (Ireland) Ltd
19 Baggot Street Lower,
Dublin 2, D02 X658, Ireland
E-mail: Ireland@arazygroup.com

Manufacturer
Blue Neem Medical Devices Pvt Ltd
Plot Nos 270 & 271, Road No 5, Harohalli Industrial Area
II Phase, Kanakapura Taluk, Ramanagara
Karnataka - 562112, India.
Phone: +91 80 29761336
E-Mail: contact@blueneem.com
Website: www.blueneem.com

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