Re-Entry Malecot Catheter

Device Description

Re-Entry Malecot Catheter is flexible tube made of polyurethane, intended to be inserted into the renal pelvis percutaneously with self-retaining Malecot winged tip to perform external drainage of urine from kidney and long tapered end inserted into ureter to perform internal drainage to bladder. Device packaging is provided with two plastic/metal retention stiffeners & urine bag connector.

For product specifications refer to the label.

Re-Entry Malecot Catheters package includes (package components may vary):

- Radiopaque catheter (with or without Hydrophilic coating)
- Plastic/metal cannulated stiffeners
- Urine bag connector

Size (Fr)	Length (cm)
8/4, 10/6, 12/6, 14/6, 16/8 & 18/8	45

Intended Purpose & Duration

Device is intended to be inserted percutaneously into the renal pelvis and to gain access to the upper urinary tract & the ureter to perform external drainage of urine from the renal pelvis and internal drainage into the bladder.

Indications For Use & Contraindications

Indicated when diagnosed to place drainage catheter percutaneously in kidney for renal drainage. No known contraindications.

Device should not be used in the conditions of Bleeding Diathesis and uncontrolled hypertension, Anticoagulant use.

Potential Complications

Potential complications include catheter occlusion and/or dislodgement, hemorrhage or hematoma, infection and pain.

Warnings And Precautions

- If a catheter has become malpositioned or if drainage ceases, the catheter should be promptly exchanged or removed.
- Do not use if product is damaged or expired and do not reuse the device.
- Use only accessories supplied along with package and do not resterilize.
- Manipulation of the product needs ultrasound or fluoroscopic guidance.
- If the device is hydrophilic coated, activate the coat by using sterile water or saline or other isotonic solution to make the catheter insertion easier.
- Dispose off the device and its components, if any, as per institutional and local regulatory guidelines for bio-hazardous medical waste disposal.
- This device is intended for use by physicians trained and experienced in diagnostic and interventional techniques, standard techniques for placement of percutaneous drainage catheters should be employed.
- Periodic evaluation of the catheter is advised; catheter should not remain indwelled for more than 30 days after deployment.

Instructions For Use

For **Catheters supplied with Hydrophilic coating**: Activate the hydrophilic coating, by wetting the catheter with sterile water or saline. For best results, keep catheter surface wet during placement.

Catheter Placement

- Perform the standard techniques for establishing a nephrostomy tract utilizing fluoroscopic control by leaving a 150 cm, 0.035 Inch guidewire in the collecting system.
- Confirm via fluoroscopy that the guidewire wire travels through all the length of the ureter and is well into the bladder.
- Dilate the musculofascial tract by progressing dilators sequentially from smallest size to the largest size to facilitate the passage of the catheter.
- Insert the Straightener provided along with the catheter and gently twist the hub to lock and close the wings and straighten the catheter.
- Pass the thin end of the catheter over the external end of the guidewire and gradually advance the catheter wings well into the collecting system and place the thin end of the catheter into the ureter.
- Confirm the position fluoroscopically then release the Straightener from the catheter to open Malecot flower in renal pelvis. Carefully remove the Straightener.
- Remove the guidewire holding the catheter in place and secure the catheter to the skin or to the fixing disc.
- Connect the catheter to the drainage bag using the urine bag connector.

Catheter Removal Procedure

- Detach the catheter from catheter fixing disc.
- Gently pull the catheter carefully to remove.
- If necessary, place a guidewire, into the kidney/collecting system to maintain access.

Product Recommendations

For catheter loading use 0.035inch wire guide only

Supply and Storage

The device is packed in primary pack made of polyethylene and the secondary package made of polyethylene and Tyvek. Supplied sterilized by Ethylene Oxide gas in peel open package. Device shall be stored at temperature range from 20°C to 40°C and away from direct sunlight.

Target Population Group

In patients who have, or might have, a urinary obstruction due to various reasons, an external drainage tube called a 'catheter' is commonly placed in the kidney in order to prevent or to temporarily relieve the obstruction

Report Any Serious Incident

If problems occur using this device or any serious incidents or adverse events, please contact our customer relations department or call to the contact provided in our manufacturing address and competent authority of the Member State in which the user and/or patient is established.

Symbol Glossory And Definations

Standard Reference:

ISO 15223-1 Medical devices — Symbols to be used with information to be supplied by the manufacturer

Symbol	Symbol Title
	Do not re-sterilize
	Do not use if package is damaged and consult instructions for use
	Consult instructions for use
	Temperature limit
	Keep away from sunlight
	Do not re-use
	Not made with natural rubber latex
	Sterilized using ethylene oxide
	Non-pyrogenic
	Keep dry
	Medical device
	Single sterile barrier system with protective packaging inside
	Does not contain DEHP
	Unique device identifier
	Date of Manufacture
	Use-by date
	Batch code

	Catalogue number
	Authorized Representative in the European Community/ European Union
	Manufacturer
	Caution
	Prescription

- This device is to be sold by/on the directions of a physician. Sterile (ETO) if package is unopened or undamaged.
- Read all the instructions before use.
- On completion of the procedure, dispose off device and its components if any as per institutional and local regulatory guidelines for bio-hazardous medical waste dispose.

Authorized Representative in the European Community
 Arazy Group (Ireland) Ltd
 19 Baggot Street Lower,
 Dublin 2, D02 X658, Ireland
 E-mail: Ireland@arazygroup.com

Manufacturer
Blue Neem Medical Devices Pvt Ltd
 Plot Nos 270 & 271, Road No 5, Harohalli Industrial Area
 II Phase, Kanakapura Taluk, Ramanagara
 Karnataka - 562112, India.
 Phone: +91 80 29761336
 E-Mail: contact@blueneem.com
 Website: www.blueneem.com

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