Suprapubic Catheter

Device Description

Suprapubic catheter is a flexible drainage tube made of polyurethane with winged petals which enables retention.

Device package includes:

- Radiopaque catheter with winged petals,
- Metal cannula with metal trocar,
- Scalpel blade and
- Urine bag connector.

Size (Fr)

8, 10, 12, 14, 16, 18 & 20

Length (cm)

25 & 30

For product specifications refer to the label.

Intended Purpose & Duration

Suprapubic catheter is intended to be inserted percutaneously into the urinary bladder and perform external drainage of urine from the urinary bladder for a period of up to 30 days.

Indications For Use & Contraindications

Indicated when diagnosed to place drainage catheter percutaneously in bladder for external drainage. Device should not be used in the conditions where in patients with previous lower abdominal surgery or Patients with acute pelvic trauma and Absence of an easily palpable or ultrasonographically localized distended urinary bladder.

Potential Complications

Potential complications include but not limited to Urinary tract infection, Bladder perforation.

Warnings And Precautions

- Manipulation of device requires Ultrasound / Fluoroscopy / Other image guiding techniques.
- Do not use if product is damaged or expired.
- This device is intended for use by physicians trained and experienced in diagnostic and percutaneous techniques, standard techniques for placement of percutaneous drainage catheters should be employed. Use aseptic techniques for deployment.
- Periodic evaluation of the device is required. If drainage ceases, the catheter should be promptly exchanged or removed.
- Use only accessories supplied along with package and do not re-sterilize.
- On completion of the procedure, dispose off the device and its components, if any, as per institutional and local regulatory guidelines for bio-hazardous medical waste disposal.
- Catheter should be irrigated on routine basis to ensure function.
- Do not reuse the device as it leads to potential complications.
- Intended for single use only.

Catheter Placement

- Identify the puncture site appropriately and anesthetize the area.
- Advance the metal cannula into the catheter tip and secure with the lure lock to close the Malecot wings, then advance the Needle through the metal cannula and secure to the hub.
- If the bladder is not already distended by urine, fill the bladder with sterile saline.
- Make an incision at the puncture site with the scalpel.
- Insert the catheter into the bladder with the help of trocar needle. Remove the trocar to observe the urine flow.
- Aspirate if necessary to confirm.
- Make sure that the catheter is stabilized, remove the needle and metal cannula carefully and then release the luer lock and rotate the hub counter clockwise to open the malecot wings.
- Secure the catheter to the catheter fixing disc or to the skin.
- Attach the urine bag connector to the catheter and allow the bladder to empty.

Catheter Removal

- Detach the catheter from catheter fixing disc or skin.
- Gently pull the catheter carefully to remove out.
- If necessary, place a guidewire, into the kidney/collecting system to maintain access.

Target Population Group

In patients who have, or might have, a urinary obstruction due to various reasons, an external drainage tube called a ‘catheter’ is commonly placed in the bladder in order to prevent or to temporarily relieve the obstruction.

Supply and Storage

The device is packed in primary pack made of polyethylene and the secondary package made of polyethylene and Tyvek. Supplied sterilized by Ethylene Oxide gas in peel open package. Device shall be stored at temperature range from 20°C to 40°C and away from direct sunlight.

Report Any Serious Incident

If problems occur using this device or any serious incidents or adverse events, please contact our customer relations department or call to the contact provided in our manufacturing address and competent authority of the Member State in which the user and/or patient is established.

Symbol Glossory And Definations

Standard Reference

ISO 15223-1 Medical devices Symbols to be used with information to be supplied by the manufacture

Symbol	Symbol Title
	Do not re-sterilize
	Do not use if package is damaged and consult instructions for use
	Consult instructions for use
	Temperature limit
	Keep away from sunlight
	Do not re-use
	Not made with natural rubber latex
	Sterilized using ethylene oxide
	Non-pyrogenic
	Keep dry
	Medical device
	Single sterile barrier system with protective packaging inside
	Does not contain DEHP
	Unique device identifier
	Date of Manufacture
	Use-by date
	Batch code

	Catalogue number
	Authorized Representative in the European Community/ European Union
	Manufacturer
	Caution
	Prescription

- This device is to be sold by/on the directions of a physician. Sterile (ETO) if package is unopened or undamaged.
- Read all the instructions before use.
- On completion of the procedure, dispose off device and its components if any as per institutional and local regulatory guidelines for bio-hazardous medical waste dispose.

Authorized Representative in the European Community
Arazy Group (Ireland) Ltd
19 Baggot Street Lower, Dublin 2,
D02 X658, Ireland
E-mail: Ireland@arazygroup.com

Manufacturer
Blue Neem Medical Devices Pvt Ltd
Plot Nos 270 & 271, Road No 5, Harohalli Industrial Area
II Phase, Kanakapura Taluk, Ramanagara
Karnataka - 562112, India.
Phone: +91 80 29761336
E-Mail: contact@blueneem.com
Website: www.blueneem.com

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