Suprapubic Drainage Catheter Set

Device Description

Suprapubic drainage catheterization inserted into the bladder percutaneously and to gain access into the urinary bladder. For product specifications refer to the label. The set contains

- Silicone foley catheter
- Plastic trocar with metallic needle tip
- Plastic peel away sheath and
- Scalpel blade.

Size (Fr)

10, 12, 14, 16 & 18

Intended Purpose

Suprapubic Drainage Catheter Set is intended to create suprapubic access to the bladder and place a catheter for bladder drainage.

Indications For Use & Contraindications

Indicated when diagnosed to place drainage catheter percutaneously in bladder for external drainage. Device should not be used in the conditions where in patients with previous lower abdominal surgery or Patients with acute pelvic trauma and Absence of an easily palpable or ultrasonographically localized distended urinary bladder.

Potential Complications

Potential complications include but not limited to Urinary tract infection, Bladder perforation.

Warnings And Precautions

- Manipulation of device requires Ultrasound / Fluoroscopy / Other image guiding techniques.
- Do not use if product is damaged or expired.
- This device is intended for use by physicians trained and experienced in diagnostic and percutaneous techniques, standard techniques for placement of percutaneous drainage catheters should be employed. Use aseptic techniques for deployment.
- If drainage ceases, the catheter should be promptly exchanged or removed.
- Use only accessories supplied along with package and do not resterilize.
- On completion of the procedure, dispose off the device and its components, if any, as per institutional and local regulatory guidelines for bio-hazardous medical waste disposal.
- Periodic evaluation of the catheter is advised.
- Catheter should be irrigated on routine basis to ensure function.
- Do not reuse the device as it leads to potential complications.
- Not intended for permanent use.

Instructions For Use

Catheter placement - Plastic trocar with sheath

- Identify the puncture site appropriately and anesthetize the area.
- If the bladder is not distended by urine, fill the bladder with sterile saline.
- Make a incision at the puncture site with the scalpel.
- Insert the plastic trocar along with peel away sheath into the bladder and observe for urine flow to confirm the access.
- Remove the plastic trocar leaving peel away sheath in place.
- Through the peel away sheath, insert the drainage catheter into the bladder and inflate the balloon present at the catheter tip.
- After inflation peel away the sheath using the holder available and remove leaving the catheter in the bladder.
- Secure the catheter to the catheter-fixing disc or to the skin.
- Attach the urine bag connector to the catheter and allow the bladder to empty.

Catheter placement - Metallic split trocar

- Identify the puncture site appropriately and anesthetize the area.
- If the bladder is not distended by urine, fill the bladder with sterile saline.
- Make an incision at the puncture site with the scalpel.
- Puncture at the site by metallic split trocar into the bladder and observe for urine flow to confirm the access.
- Through the split trocar, insert drainage catheter into the bladder. Then split the trocar half away by holding on both the sides of hub of trocar. Then remove the parted trocar slits and inflate the balloon present inside the bladder by sterile water or saline. Do not use air, any gas to inflate the balloon.
- Secure the catheter to skin by catheter-fixing disc or surgical tape.
- Attach the urine bag connector to the catheter and allow the bladder to empty.

Catheter Removal

- Detach the catheter from catheter fixing disc or skin.
- Deflate the balloon of the catheter.
- Gently pull the catheter carefully to remove.

Target Population Group

In patients who have, or might have, a urinary obstruction due to various reasons, an external drainage tube called a 'catheter' is commonly placed in the bladder in order to prevent or to temporarily relieve the obstruction

Supply and Storage

The device is packed in blister tray sealed with medical grade Tyvek pouch. Supplied sterilized by Ethylene Oxide gas in peel open package. Device shall be stored at temperature range from 20OC to 40OC and away from direct sunlight.

Report Any Serious Incident

If problems occur using this device or any serious incidents or adverse events, please contact our customer relations department or call to the contact provided in our manufacturing address and competent authority of the Member State in which the user and/or patient is established.

Symbol Glossary And Definitions

Standard Reference:

ISO 15223-1 Medical devices - Symbols to be used with information to be supplied by the manufacturer

Symbol	Symbol Title
	Do not re-sterilize
	Do not use if package is damaged and consult instructions for use
	Consult instructions for use
	Temperature limit
	Keep away from sunlight
	Do not re-use
	Not made with natural rubber latex
	Sterilized using ethylene oxide
	Non-pyrogenic
	Keep dry
	Medical device
	Single sterile barrier system
	Does not contain DEHP
	Unique device identifier
	Date of Manufacture
	Use-by date
	Batch code

	Catalogue number
	Authorized Representative in the European Community/ European Union
	Manufacturer
	Caution
	Prescription

- This device is to be sold by/on the directions of a physician. Sterile (ETO) if package is unopened or undamaged.
- Read all the instructions before use.
- On completion of the procedure, dispose off device and its components if any as per institutional and local regulatory guidelines for bio-hazardous medical waste dispose.

Authorized Representative in the European Community
 Arazy Group (Ireland) Ltd
 19 Baggot Street Lower, Dublin
 2, D02 X658, Ireland
 E-mail: Ireland@arazygroup.com

Manufacturer
Blue Neem Medical Devices Pvt Ltd
 Plot Nos 270 & 271, Road No 5, Harohalli Industrial Area
 II Phase, Kanakapura Taluk, Ramanagara
 Karnataka - 562112, India.
 Phone: +91 80 29761336
 E-Mail: contact@blueneem.com
 Website: www.blueneem.com

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