Ureteral Balloon Dilator

Device Description

Balloon Dilator contains a distal balloon which can show inflation and deflation property. Balloon have metal radiopaque bands on proximal and distal ends of balloon. Device package contains (contents may vary). Refer to product label for specifications.

- Balloon dilator with syringe
- Balloon dilator with inflation device

Intended Purpose

Intended for dilation of narrow ureteral segments through natural orifice prior to a diagnostic or interventional procedure.

Indications For Use & Contraindications

Do not use this device to treat obstruction that is due to extrinsic forces. This device is contraindicated in the presence of conditions which create unacceptable risks during catheterization.

Specifications

Balloon expansion diameter: 05 mm
 Balloon length: 04 cm
 Nominal pressure: 07 atm
 Rate burst pressure: 14 atm
 Use guide wire: 0.035 inch

Potential Complications

Potential complications Include but not limited to Infection, Tissue perforation, Tissue trauma, Hematuria.

Warnings And Precautions

- This device is intended for use by physicians trained and experienced in techniques for balloon dilation or catheterization and this device should not be used without comprehensive knowledge of the indications, techniques and risks of the procedure.
- Always inflate the balloon with sterile liquid. Do not inflate the balloon with any kind of air or carbon dioxide or other gas. Do not over inflate or exceed maximum rated burst pressure, as excessive pressure may cause the balloon to rupture.
- Refer to product label for appropriate inflation details.
- If loss of pressure within the balloon occurs during inflation or if balloon ruptures during dilatation, immediately discontinue the procedure. Deflate the balloon carefully. Do not attempt to reinflate the balloon. Carefully withdraw the device if failure found during procedure.
- Do not exceed the maximum rated burst pressure (listed on label) for this balloon device. Refer label for device specifications.
- Use only compatible wire guide as specified on the device label.
- Do not use if product is damaged or expired, reuse and resterilization of this device leads to loss in structural integrity and performance.
- On completion of the procedure, dispose off the device as per local regulatory guidelines for bio-hazardous medical waste disposal

Instruction of Use

- Remove the sleeve from the device and discard.
- Lock the syringe on to the side arm stop cock of the balloon catheter and aspirate until all air is removed and close the stop cock.
- Using imaging guidance, pass a guidewire of 0.035 Inch to the desired distance of the tract to made.
- Pass the balloon dilator over the guidewire and position the balloon with the help of radiopaque markers under imaging guidance, Please Ensure proper placement of the balloon at the site of dilation.
- Fill the syringe or use the inflation device if provided, with contrast medium and eliminate air from the syringe or inflation device.
- Open the stop cock on the dilator, connect the syringe or inflation device and inflate the balloon with the contrast media.

Target Population Group

Adult Male and Female population groups.

CAUTION: Inflate the balloon only with sterile medium. Do not inflate with air or any gas. To ensure proper regulation of balloon pressure, use of balloon inflation device with pressure gage is recommended. Do not exceed the recommended max pressure mentioned on the product label.

- On reaching desired balloon pressure, close the stop cock and leave the inflated balloon at place for at least 30 seconds to allow adequate dilation.
- Once desired dilation is achieved, open the stop cock and aspirate the liquid from the balloon using the syringe or inflation device.
- For removal of the balloon, aspirate the balloon completely and gently rotate the dilator shaft counter- clockwise during withdrawal.

Supply and Storage

The device is assembled in protective dispenser tube and the secondary package made of polyethylene and Tyvek. Supplied sterilized by Ethylene Oxide gas in peel open package. Device shall be stored at temperature range from 20°C to 40°C and away from direct sunlight.

Report Any Serious Incident

If problems occur using this device or any serious incidents or adverse events, please contact our customer relations department or call to the contact provided in our manufacturing address and competent authority of the Member State in which the user and/or patient is established.

Symbol Glossary And Definitions

Standard Reference:

ISO 15223-1 Medical devices - Symbols to be used with information to be supplied by the manufacturer

Symbol	Symbol Title
	Do not re-sterilize
	Do not use if package is damaged and consult instructions for use
	Consult instructions for use
	Temperature limit
	Keep away from sunlight
	Do not re-use
	Not made with natural rubber latex
	Sterilized using ethylene oxide
	Non-pyrogenic
	Keep dry
	Medical device
	Single sterile barrier system with protective packaging inside
	Does not contain DEHP
	Unique device identifier
	Date of Manufacture
	Use-by date
	Batch code

	Catalogue number
	Authorized Representative in the European Community/ European Union
	Manufacturer
	Caution
	Prescription

- This device is to be sold by/on the directions of a physician. Sterile (ETO) if package is unopened or undamaged.
- Read all the instructions before use.
- On completion of the procedure, dispose off device and its components if any as per institutional and local regulatory guidelines for bio-hazardous medical waste dispose.

Authorized Representative in the European Community
Arazy Group (Ireland) Ltd
19 Baggot Street Lower, Dublin 2,
D02 X658, Ireland
E-mail: Ireland@arazygroup.com

Manufacturer
Blue Neem Medical Devices Pvt Ltd
Plot Nos 270 & 271, Road No 5, Harohalli Industrial Area
II Phase, Kanakapura Taluk, Ramanagara
Karnataka - 562112, India.
Phone: +91 80 29761336
E-Mail: contact@blueneem.com
Website: www.blueneem.com

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Auktoriserad representant i Europeiska gemenskapen Arazy Group (Ireland) Ltd