Ureteral Catheter

Device Description

Ureteral Catheter is a flexible tube-like device made from PEBA for introduction into urinary tract through an endoscope. Device is available in variety of French sizes, lengths and distal tip configurations. Choice of device is based on clinical situation. Refer to label for technical specifications.

Size (Fr)	Length (cm)	Model
3, 4, 5, 6, 7 & 8	70	- Straight Open End & Closed End - Angled Tip Open End & Closed End - Cone Tip Open End & Closed End

Intended Purpose & Duration

Ureteral Catheters are used for catheterization of the urinary tract during interoperative procedures, includes the following operations:

- Injection of contrast media such as to perform retrograde pyelogram.
- Navigation of a tortuous ureter.
- Access, drainage of fluids and advancement or exchange of ancillary devices including guidewires (for open end catheters only) for a period of up to 06 hours.

Indications For Use & Contraindications

Contraindicated for use in patients unsuitable for catheterization

Potential Complications

Potential complications Include but not limited to haemorrhage, Edema, perforation of the kidney, ureter avulsion renal pelvis, ureter and/or bladder, urinary tract infection.

Warnings And Precautions

- Manipulation of device requires Ultrasound / Fluoroscopy / Other image guiding techniques.
- Do not use if product is damaged or expired, reuse and resterilization of this device leads to potential complications and loss in structural integrity.
- This device is intended for use by physicians trained and experienced in diagnostic and endoscopic techniques, standard techniques for placement of endoscopic devices should be employed.
- On completion of the procedure, dispose off the device and its components, if any, as per institutional and local regulatory guidelines for bio-hazardous medical waste disposal.
- Not intended for permanent use.
- Use only compatible endoscopic devices.
- Avoid bending or kinking of catheter prior to placement. Doing so could damage the integrity of the catheter and result in patient injury.

Catheter Placement Procedure

Instructions For Use

For Closed End - Supplied with Stylet (Placed without the guidewire)

1. Place the endoscope into bladder and view the ureteral orifice, Choose the appropriate catheter type based on patient anatomy.
2. The catheter adapter is secured to the ureteral catheter.
3. Ensure that the stylet is inserted into the catheter.
4. Under fluoroscopic vision the catheter tip is positioned within the ureteral orifice through the endoscope.
5. Remove the stylet carefully.
6. Contrast media can then be injected through the ureteral catheter into the ureter.

For Open End (Placed with the guidewire)

1. Place the endoscope into the bladder and view the ureteral orifice. Choose the appropriate catheter type based on patient anatomy.
2. Under fluoroscopic guidance, place the guidewire within the ureteral orifice.
3. The ureteral catheter is inserted over the guidewire and advanced into the ureter.
4. Remove the guidewire carefully.
5. Secure the adapter to the catheter.
6. Contrast media can then be injected through the ureteral catheter into the ureter.

NOTE: This Device is detectable by X-Ray or other means of fluoroscopy.

Catheter Removal Procedure

Once the ancillary device is placed or in position or contrast media is injected, detach the catheter hub manually and then gently withdraw the catheter.

Product Recommendations

Recommended wire guide for the usage of this product is as follows:

- 03Fr : 0.018inch
- 04Fr : 0.025inch
- 05 Fr : 0.035inch
- 06, 07 & 08Fr : 0.038 inch

Supply and Storage

The device is packed in primary pack made of polyethylene and the secondary package made of polyethylene and Tyvek. Supplied sterilized by Ethylene Oxide gas in peel open package. Device shall be stored at temperature range from 20°C to 40°C and away from direct sunlight.

Report Any Serious Incident

If problems occur using this device or any serious incidents or adverse events, please contact our customer relations department or call to the contact provided in our manufacturing address and competent authority of the Member State in which the user and/or patient is established.

Target Population Group

Intended for both adult and pediatric patients for accessing and catheterizing the urinary tract allowing for procedures like contrast medium delivery, fluid drainage, and irrigation. Catheter choice is based on physicians preference and clinical situation.

Symbol Glossory And Defanitions

Standard Reference:

ISO 15223-1 Medical devices — Symbols to be used with information to be supplied by the manufacturer

Symbol	Symbol Title
	Do not re-sterilize
	Do not use if package is damaged and consult instructions for use
	Consult instructions for use
	Temperature limit
	Keep away from sunlight
	Do not re-use
	Not made with natural rubber latex
	Sterilized using ethylene oxide
	Non-pyrogenic
	Keep dry
	Medical device
	Single sterile barrier system with protective packaging inside
	Does not contain DEHP
	Unique device identifier
	Date of Manufacture
	Use-by date
	Batch code

	Catalogue number
	Authorized Representative in the European Community/ European Union
	Manufacturer
	Caution
	Prescription

- This device is to be sold by/on the directions of a physician. Sterile (ETO) if package is unopened or undamaged.
- Read all the instructions before use.
- On completion of the procedure, dispose off device and its components if any as per institutional and local regulatory guidelines for bio-hazardous medical waste dispose.

Authorized Representative in the European Community
Arazy Group (Ireland) Ltd
19 Baggot Street Lower,
Dublin 2, D02 X658, Ireland
E-mail: Ireland@arazygroup.com

Manufacturer
Blue Neem Medical Devices Pvt Ltd
Plot Nos 270 & 271, Road No 5, Harohalli Industrial Area
II Phase, Kanakapura Taluk, Ramanagara
Karnataka - 562112, India.
Phone: +91 80 29761336
E-Mail: contact@blueneem.com
Website: www.blueneem.com

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