en Instructions for Ureteral Dilators

Device Description

Ureteral Dilators are long tube-like device made of medical grade plastic material, with a tapered tip at the distal end. The Dilators are radiopaque for imaging guidance. Device is available in individual dilators and dilator sets. Refer to the product label for specification.

Ureteral Dilator		Ureteral Dilator Set
Size (Fr)	Length (cm)	6-16 Fr (6, 8, 10, 12, 14 & 16 Fr)
6, 8, 10, 12, 14 & 16	70	

Intended Purpose

Intended used for dilation of the ureter through natural orifice prior to a diagnostic or interventional procedure.

Indications For Use & Contraindications

None known.

Potential Complications

Potential complications Include but not limited to Edema, Irritation, inflammation, Ureteral avulsion.

Warnings And Precautions

- Use with caution when advancing the dilators to avoid tissues and orifice damages. Use only compatible wireguide.
- Do not use if product is damaged or expired, reuse and resterilization of this device leads to loss in structural integrity and performance.
- This device is intended for use by trained physicians and this device should not be used without comprehensive knowledge of the indications, techniques and risks of the procedure.
- On completion of the procedure, dispose off the device as per local regulatory guidelines for bio-hazardous medical waste disposal.
- Imaging guidance is required for use of this device.

Instruction of Use

- Before use, immerse the dilator in sterile water or isotonic solution to allow the hydrophilic surface to absorb water and become lubricious makes its passage easier under normal conditions (for hydrophilic coated dilators only) (refer to device label).
- Under x-ray or fluoroscopic control, endoscopically pass the guidewire to the desired distance into the ureter under imaging guidance or through endoscope.
- Progressing from the smallest to the largest appropriate size, pass the radiopaque dilators over the guidewire while maintaining the Guidewire's position.

Target Population Group

Adult male and female individuals

Supply and Storage

The device is packed in primary pack made of polyethylene and the secondary package made of polyethylene and Tyvek. Supplied sterilized by Ethylene Oxide gas in peel open package. Device shall be stored at temperature range from 20°C to 40°C and away from direct sunlight.

Report Any Serious Incident

If problems occur using this device or any serious incidents or adverse events, please contact our customer relations department or call to the contact provided in our manufacturing address and competent authority of the Member State in which the user and/or patient is established.

Symbol Glossory And Definations

Standard Reference:

ISO 15223-1 Medical devices — Symbols to be used with information to be supplied by the manufacturer

Symbol	Symbol Title
	Do not re-sterilize
	Do not use if package is damaged and consult instructions for use
	Consult instructions for use
	Temperature limit
	Keep away from sunlight
	Do not re-use
	Not made with natural rubber latex
	Sterilized using ethylene oxide
	Non-pyrogenic
	Keep dry
	Medical device
	Single sterile barrier system with protective packaging inside
	Does not contain DEHP
	Unique device identifier
	Date of Manufacture
	Use-by date
	Batch code

	Catalogue number
	Authorized Representative in the European Community/ European Union
	Manufacturer
	Caution
	Prescription

- This device is to be sold by/on the directions of a physician. Sterile (ETO) if package is unopened or undamaged.
- Read all the instructions before use.
- On completion of the procedure, dispose off device and its components if any as per institutional and local regulatory guidelines for bio-hazardous medical waste dispose.

Authorized Representative in the European Community
Arazy Group (Ireland) Ltd
19 Baggot Street Lower,
Dublin 2, D02 X658, Ireland
E-mail: Ireland@arazygroup.com

Manufacturer
Blue Neem Medical Devices Pvt Ltd
Plot Nos 270 & 271, Road No 5, Harohalli Industrial Area
II Phase, Kanakapura Taluk, Ramanagara
Karnataka - 562112, India.
Phone: +91 80 29761336
E-Mail: contact@blueneem.com
Website: www.blueneem.com

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