

Amplatz Ultra Stiff[®] Guidewire



Types

- J Tip

Size (in inch)

0.032

Length (in cm)

80

Intended Use:

Used to gain access to desired anatomical location and to assist in placement of medical devices during diagnostic or interventional procedures (PCNL).

Key Features:

- PTFE coated stainless steel mandril with 5 cm atraumatic floppy tip
- Smooth transition from stiff mandril to the flexible tip
- The stiffness of the mandril design provides support for the advancement of devices such as dilators, catheters and instruments