

Grasping Forceps (Triprong)



Intended Use:

Used to remove the stones and retrieval of foreign objects under direct vision.

URS	
Size (in Fr)	2.5, 3, 4 & 5
Length (in cm)	60

PCNL	
Size (in Fr)	7, 8, 9 & 12
Length (in cm)	40

Cystoscope	
Size (in Fr)	2.5, 3, 4 & 5
Length (in cm)	60

Flexible Graspers	
Size (in Fr)	2.5, 3, 4 & 5
Length (in cm)	40, 60, 70 110, 115 & 120