


## Hydro Twister<sup>®</sup> Guidewire

### Intended Use:

Used in placement of medical devices while gaining easy access to desired anatomical location



Type	Size (in inch)	Length (in cm)
Straight Tip	0.025, 0.032, 0.035 & 0.038	150

### Key Features and Benefits:

- **Hydrophilic coated tip and sleeve:** Easy, low-friction advancement, excellent visualization, and maneuverability during procedures (URS and RIRS)
- **Atraumatic floppy tip:** Facilitates easy insertion and access
- **Kink-resistant nitinol core:** Reliable deployment of devices
- **Single wire compatibility:** Avoid exchange of wires; Less procedural time; Easy access in tortuous anatomy; Enable access and stenting with a single wire