


Hydro Twister[®] Guidewire



Type
Straight Tip

Size (in inch)
0.025, 0.032, 0.035
& 0.038

Length (in cm)
150

Intended Use:

Used in placement of medical devices while gaining easy access to desired anatomical location. (URS & RIRS)

Key Features:

- Hydrophilic coated tip and blue & white sleeve provides
 - Easy, low friction advancement with excellent visualization
 - Maneuverability during diagnostic and interventional procedures
- Atraumatic floppy tip facilitates easy insertion and access
- Is constructed with kink resistant Nitinol core with appropriate stiffness for deployment of devices

Advantages:

- Avoid exchange of wires
- Less procedural time
- Easy access in tortuous anatomy
- Access and stenting both can be performed with a single wire