


## Hydro Twister<sup>®</sup> Plus Guidewire

### Intended Use:

Used to gain access to desired anatomical location and to assist in placement of medical devices



Type	Size (in inch)	Length (in cm)
Straight Tip	0.025, 0.032, 0.035 & 0.038	150

### Key Features and Benefits:

- **Hydrophilic coated tip:** Easy, low-friction advancement and excellent visualization
- **Enhanced maneuverability:** Facilitates diagnostic and interventional procedures
- **Atraumatic floppy tip:** Easy insertion and access
- **Kink-resistant nitinol core:** Appropriate stiffness for device deployment
- **Double-sided floppy tip:** Prevents scope damage during multiple insertions
- **Single wire access:** Avoids exchange of wires; access and stenting can both be performed
- **Easy access in tortuous anatomy:** Less procedural time