

Hydro Twister[®] Plus



Type
Straight Tip

Size (in inch)
0.025, 0.032, 0.035
& 0.038

Length (in cm)
150

Intended Use:

Used to gain access to desired anatomical location and to assist in placement of medical devices. (URS & RIRS)

Key Features:

- Hydrophilic coated tip and blue & white sleeve provides
 - Easy, low friction advancement with excellent visualization
 - Maneuverability during diagnostic and interventional procedures
- Atraumatic floppy tip facilitates easy insertion and access
- Is constructed with kink resistant Nitinol core with appropriate stiffness for deployment of devices
- Double side floppy tip prevents the damage of scopes due to multiple insertions during the procedure

Advantages:

- Avoid exchange of wires
- Less procedural time
- Easy access in tortuous anatomy
- Access and stenting both can be performed with a single wire