

HYDRO TWISTER[®] GUIDEWIRE

Used in placement of medical devices while gaining easy access to desired anatomical location (URS & RIRS).

Key Features

- Hydrophilic coated tip and blue & white sleeve provides
 - Easy, low friction advancement with excellent visualization
 - Maneuverability during diagnostic and interventional procedures
- Atraumatic floppy tip facilitates easy insertion and access
- Is constructed with kink resistant Nitinol core with appropriate stiffness for deployment of devices

Advantages

- Avoid exchange of wires
- Less procedural time
- Easy access in tortuous anatomy
- Serves two purpose - access and stenting with single wire

Type

Straight Tip

Size (in inch)

0.025, 0.032, 0.035 & 0.038

Length (in cm)

150

