

HYDRO TWISTER[®] PLUS

Used to gain access to desired anatomical location and to assist in placement of medical devices (URS & RIRS).

Key Features

- Is constructed with kink resistant nitinol core with appropriate stiffness for deployment of devices such as stents during a procedure
- Atraumatic floppy tip of 5 cm facilitates easy insertion and access
- Hydrophilic coated tip and blue & white sleeve provides easy, low friction advancement with excellent visualization and maneuverability during diagnostic and interventional procedures
- Double side floppy tip prevents the damage of scopes due to multiple insertions during the procedures

Advantages

- Avoid exchange of wires
- Less procedural time
- Easy access in tortuous anatomy
- Serves two purpose - access and stenting with single wire

Type

Straight Tip

Size (in inch)

0.025, 0.032, 0.035 & 0.038

Length (in cm)

150

