

Initial Puncture Needle Two Part



Type
Two Part

Size (in G)
14, 16 & 18

Length (in cm)
10, 15, 20, 23 & 25

Intended Use:

Intended to make the initial percutaneous puncture to access the kidney and to enable the introduction of a guidewire during urological procedures (PCN, PCNL, Mini PERC & Micro PERC).

Key Features:

- Siliconised needle for less trauma
- Echogenic tip for better visualization during ultrasound
- 1 cm marking on the shaft

Advantage of Locking hub IP Needle:

During the procedure locking hub helps in holding the inner needle with outer cannula together, which helps in easy initial puncture