

INITIAL PUNCTURE NEEDLE - TWO PART

Intended to make the initial percutaneous puncture to access the kidney and to enable the introduction of a guidewire during urological procedures (PCN, PCNL, Mini PERC & Micro PERC).

Key Features

- Siliconised needle for less trauma
- Echogenic tip for better visualization during ultrasound
- 1 cm marking on the shaft

Advantage of Locking Hub IP Needle

During the procedure the locking hub helps in holding the inner needle with outer cannula together, which helps in easy initial puncture

Type

Two Part

Size (in G)

14, 16 & 18

Length (in cm)

10, 15, 20, 23 & 25

