

Initial Puncture Needle - Two Part

Intended Use:

Intended to make the initial percutaneous puncture to create access and to enable the introduction of guidewire during percutaneous interventional procedures



Type	Size (in G)	Length (in cm)
Two part	14, 16 & 18	10, 15, 20, 23 & 25

Key Features and Benefits:

- **Siliconized needle:** Less trauma during the initial puncture
- **Echogenic tip:** Better visualization during ultrasound improving accuracy and safety
- **1cm shaft marking:** Helps in easy initial puncture
- **Locking hub:** Helps in holding the inner needle with outer cannula together, which helps in easy initial puncture