

PCN Catheter



Package Includes

PCN Catheter, metal or plastic straightener & urine bag connector

Size (in Fr)

6, 7, 8, 8.5, 9, 10, 12, 14
16 & 18

Length (in cm)

20, 22, 25, 30, 35 & 40

Intended Use:

PCN Catheter is intended to be inserted percutaneously into the renal pelvis to perform external drainage of urine.

Key Features:

- The pigtail tip configuration helps in the retention of the distal tip in the renal pelvis
- Eye holes for better drainage
- Also available with Hydrophilic coating for easy insertion