


PCN Set - Pigtail (Nephrostomy Catheter Set)

Intended Use:

Used for percutaneous placement of a pigtail catheter in the renal pelvis to perform external drainage of urine.



Set Includes

- IP Needle
- Guidewire J tip PTFE coated
- Fascial Dilators
- Urine bag connector
- PCN Catheter
- Scalpel blade
- Catheter fixing disc

Size (in Fr)

6, 7, 8, 8.5, 9, 10, 12, 14, 16 & 18

Length (in cm)

20, 22, 25, 30 & 35

Benefit:

Kit provides entire devices for performing Nephrostomy procedure