

Praveza[®] Ureteral Access Sheath



Sizes available for all types of flexible scopes

Storz, Wolf & Single use flexible scopes

Size (in Fr)

ID	9.5	10	11	12
OD	11.5	12	13	14

Length (in cm)

28, 35, 45 & 55

Intended Use:

Used to establish a conduit during urological procedures facilitating the passage for Ureteroscope/ Endoscope and other surgical instruments into the urinary tract (RIRS).

Key Features:

- Kink resistant – reinforced with a stainless steel coil
- Designed to provide stability and strength through tortuous anatomy
- Tapered tip dilator with hydrophilic coating provides easy access
- Transparent tip with platinum ring helps in better visualization and placement under fluoroscopy
- Sizes available for all types of flexible scopes
- Seamless transition between dilator and sheath for easy and atraumatic insertion