

Praveza[®]

Ureteral Access Sheath

Intended Use:

Used to establish a conduit during urological procedures facilitating the passage for Ureteroscope/ Endoscope and other surgical instruments into the urinary tract (RIRS).



Sizes available for all types of flexible scopes
Storz, Wolf and Single use flexible scopes

Size (in Fr)

ID	9.5	10	11	12
OD	11.5	12	13	14

Length (in cm)
28, 35, 45 & 55

Key Features and Benefits:

- Kink resistant – reinforced with a stainless steel coil
- Designed to provide stability and strength through tortuous anatomy
- Tapered tip dilator with hydrophilic coating provides easy access
- Transparent tip with platinum ring helps in better visualization and placement under fluoroscopy
- Sizes available for all types of flexible scopes
- Seamless transition between dilator and sheath for easy and atraumatic insertion