


PRAVEZA[®] URETERAL ACCESS SHEATH



Types

- Storz
- Wolf
- Single use flexible scopes

Size (in Fr)

ID	9.5	10	11	12
OD	11.5	12	13	14

Length (in cm)

28, 35, 45 & 55

Intended Use:

Used to gain smooth and easy access to desired anatomical location and to assist in placement of medical devices.

Key Features:

- Kink resistant : reinforced with a stainless steel coil designed to provide stability and strength through tortuous anatomy
- Tapered tip dilator with hydrophilic coating provides easy access
- Transparent tip with platinum ring helps in better visualization and placement of sheath under fluoroscopy
- Sizes available for all types of flexible scopes (Storz, Wolf & Single use flexible scopes)
- Seamless transition between dilator and sheath for easy and atraumatic insertion