

## Suprapubic Catheter



### Package Includes

Suprapubic Catheter, metal straightner, metal trocar and urine bag connector

### Size (in Fr)

8, 10, 12, 14, 16, 18 & 20

### Length (in cm)

25 & 30

### Intended Use:

Suprapubic Catheter is intended to be inserted percutaneously into the urinary bladder and perform external drainage of urine from the urinary bladder.

### Key Features:

- Integral flower design matches with the catheter size when closed and cause less trauma during insertion
- Trocar and catheter designed to prevent buckling of catheter during the procedure