

# SUPRAPUBIC CATHETER

Suprapubic catheter is intended to be inserted percutaneously into the urinary bladder and perform external drainage of urine from the urinary bladder.

## Key Features

- Integral flower design matches with the catheter size when closed and cause less trauma during insertion
- Trocar and catheter designed to prevent buckling of catheter during the procedure

## Package Includes

Suprapubic Catheter, metal straightner, metal trocar and urine bag connector

## Size (in Fr)

8, 10, 12, 14, 16, 18 & 20

## Length (in cm)

25 & 30

