

# SUPRAPUBIC DRAINAGE CATHETER SET (SDS)

Suprapubic drainage catheter set is intended to create suprapubic access to the bladder and place a catheter for bladder drainage.

## Advantages

- Silicone foley catheter for patient comfort and safety
- Specially designed trocar with needle tip enhances smooth puncture
- Peel away sheath forms a conduit for catheter placement and can be peeled easily to remove from the catheter
- Seamless transition between trocar and peel away sheath for easy and atraumatic puncture

## Set Includes

Foley catheter, trocar with peel-away sheath and scalpel blade

## Size (in Fr)

10, 12, 14, 16 & 18

