

## Twister Guidewire



### Type

Straight Tip

### Size (in inch)

0.025, 0.032, 0.035  
& 0.038

### Length (in cm)

150

### Intended Use:

Used to gain access to desired anatomical location and to assist in placement of medical devices. (URS procedure, for ureteral stenting).

### Key Features:

- Is constructed with kink resistant Nitinol core with appropriate stiffness for deployment of devices
- Atraumatic floppy tip facilitates easy insertion and access
- Blue and White sleeve provides excellent visualization and easy maneuvering