

# TWISTER GUIDEWIRE

Used to gain access to desired anatomical location and to assist in placement of medical devices (URS procedure, for ureteral stenting).

## Key Features

- Is constructed with kink resistant nitinol core with appropriate stiffness for deployment of devices such as stents during a procedure
- Atraumatic floppy tip of 5 cm facilitates easy insertion and access
- Blue & White sleeve provides excellent visualization and easy maneuvering during diagnostic and interventional procedure

## Type

Straight Tip

## Size (in inch)

0.025, 0.032, 0.035 & 0.038

## Length (in cm)

150

