


## Mono J Stent



### Types

- Open End
- Closed End

### Size (in Fr)

4.7, 5, 6, 7 & 8

### Length (in cm)

70

### Also available in set, which Includes

Mono J stent & guidewire

### Intended Use:

Used for temporary drainage and irrigation of urine from the kidney (Pigtail ureteral catheter).

### Key Features:

- Pigtail loop helps retention in renal pelvis
- Holes in shaft helps in better drainage
- Used to inject contrast
- Used for short term external drainage to avoid stenting